# ForaCare® D20 No-Coding

Blood Glucose Plus Blood Pressure Monitoring System

## Dear ForaCare D20 No-Coding System Owner:

Thank you for choosing ForaCare D20 No-Coding Blood Glucose Plus Blood Pressure Monitoring System to take control of your diabetes. Being a home care device, this system is helpful in tracking your blood glucose and blood pressure values.

The Quick Start User Guide will assist you in setting up the system and testing your blood glucose levels. Please read it carefully before use. For complete information on the system's use and features, please refer to the Owner's Manual.

**Test Strip Vial**

**Test Strips**

**Lancets**


**Control Solution**

Please note: Use only FORA D20 test strips with ForaCare D20 No-Coding System meter.

**Clear Cap**  
(for obtaining blood from places other than your fingertip)

**Lancing Device**

**Lancing device, lancets and clear cap may not be included.**  
If your lancing device differs from the one shown here, please refer to the manufacturer's manual to ensure proper usage.  
Please Note: Your lancing device may differ from the one shown in here.

## Quick Start User Guide

For self-testing. Read instructions before use. For *In vitro* diagnostic use.

Distributed by ForaCare

Fora Care Inc.  
893 Patriot Dr., Suite D, Moorpark, CA 93021

Products made in Taiwan

Toll Free:

1-888-307-8188 (8:30am-5:00pm PST, Mon.-Fri.)

1-866-469-2632 (24 Hours, 7 days/week)

www.foracare.com

315-3263100-011


### Front View of Meter

- ① Pressure Cuff
- ② On/Off Button
- ③ C Button
- ④ M Button
- ⑤ Test Slot

- ⑥ DC Adapter Port
- ⑦ PC Link Port
- ⑧ Display Screen
- ⑨ S Button

### Rear View of Meter

- ⑩ Battery Compartment
- ⑪ Battery Cover

## 1. Setting the Monitor

Press **S** to start setting. Year flashes first.

Press **M** to change values. Hold down **M** to increase the number more quickly. Press **S** to move on to the next step.

## 1. Setting the Monitor


Date flashes. Follow **1.2** to make changes. The sequence of settings is:

- |                            |   |
|----------------------------|---|
| (1) date                   | (4) memory deletion                     |
| (2) time                   | (5) voice volume                        |
| (3) unit of blood pressure | (6) language (L1: English, L2: Spanish) |

The meter will automatically turn off when the setting is complete. For memory deletion, please refer to the Owner's Manual for details.


## 2. Setting the Lancing Device

Twist and pull off the cap of the lancing device.

Insert a lancet into the lancet holder and push down firmly until it is fully secured.

Twist the protective disk off the lancet.

Replace and twist the cap to the right to close the lock as shown in the picture until you hear the click. When the lock is properly closed lock indicator is in the extreme right position.

## 2. Set the Lancing Device

Holding the cap select on the depth of penetration by turning the adjustable tip in either direction so that the arrow on the cap points to the desired depth.

Pull the cocking control back until it clicks. You will see a color change inside the release button when it is ready.

► If it does not click, the device may have been cocked when the lancet was inserted.

## 3. Testing Your Blood Glucose

Tear off the security label (Tamper Seal). When opening the vial for the first time, record the date on the label.

Wash and dry your hands thoroughly. Take a test strip out and then close the vial tightly.


Insert a test strip to turn on the monitor. The sequences of LCD display are: “ → flashing “” with date and time.


Take the pre-set lancing device. Hold the lancing device firmly against the puncture site. Press the release button to obtain a drop of blood.

## 3. Testing Your Blood Glucose

Gently squeeze around the punctured area to obtain enough blood.

While “” is flashing on the display, hold and touch the absorbent hole of the test strip to the drop of blood until the confirmation window is completely filled.

After the monitor counts down to 0, your test result will appear along with the date and time and will be automatically stored in the memory.

Remove the used test strip. Dispose of the used strip and lancet carefully.

For testing with control solution, please refer to the Owner's Manual for detailed information.

## 4. Reading Your Result

Your result will display the following messages:

| Appears        | When glucose              |
|----------------|---------------------------|
| <b>Lo</b>      | < 20 mg/dL (1.1 mmol/L)   |
| <b>KETONE?</b> | ≥ 240 mg/dL (13.3 mmol/L) |
| <b>Hi</b>      | > 600 mg/dL (33.3 mmol/L) |

Reference values for blood glucose:

| Time of day             | Normal plasma glucose range for people without diabetes |
|-------------------------|---|
| Fasting and before meal | < 100 mg/dL (5.6 mmol/L)                                |
| 2 hours after meals     | < 140 mg/dL (7.8 mmol/L)                                |

Source:

American Diabetes Association (2012). Diabetes Care, 35 (Suppl 1): S64-71.

Please work with your doctor to determine a target range that works best for you.

## 5. Testing Your Blood Pressure

Sit down for at least 10 minutes before measuring.

Stretch your left arm forward with your palm facing up. Slide the cuff onto your left arm, above your elbow. The red line on the edge of the cuff should be approximately 0.8 to 1.2 inch (2 to 3 cm) above your elbow.


Arrange tubing over the main arteries inside of your arm. Press .


Relax and make sure the cuff is at the same height as your heart.

## 5. Testing Your Blood Pressure

Measurement is in progress: Start inflation → Stop inflation when cuff pressure is enough → Cuff pressure deflates steadily

The monitor will display the systolic pressure, diastolic pressure, and heart rate along with the date and time. Press to switch off.

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