
FORA[®] GD50

Blood Glucose Monitoring System

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Meter Serial Number (See label on back of meter) /

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Foracare, Inc. Online Warranty Registration and Detail Warranty Coverage at www.foracare.com/support.html

PLACE
STAMP HERE
COLOQUE
ESTAMPILLA
AQUÍ

ForaCare, Inc.

893 Patriot Dr., Suite D

Moorpark, CA 93021 USA

Dear FORA GD50 System Owner:

Thank you for purchasing the **FORA GD50** Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read the following contents thoroughly and carefully.

This system is intended to be used by healthcare professionals for multiple patients in a professional healthcare setting (i.e. hospitals, doctor's office, etc.)

If you have questions, just call Foracare Customer Care Service Center toll-free at 1-888-307-8188 (7:00 am - 6:00 pm PST, Mon. - Fri.) or visit www.foracare.com.

IMPORTANT SAFETY INSTRUCTIONS

READ BEFORE USE

- Users need to adhere to Standard Precautions when handling or using this device. All parts of the glucose monitoring system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and healthcare professionals. For more information, refer to “Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings 2007”, <http://www.cdc.gov/hicpac/2007ip/2007isolationprecautions.html>.
 - The meter should be disinfected after use on each patient. This Blood Glucose Monitoring System may only be used for testing multiple patients when Standard Precautions and the manufacturer's disinfection procedures are followed.
 - Only auto-disabling, single use lancing devices may be used with this device.
1. Use this device ONLY for the intended use described.
 2. Do NOT use accessories which are not specified by the manufacturer.
 3. Do NOT use the device if it is not working properly or if it is damaged.
 4. This device does NOT serve as a cure for any symptoms or diseases.
The data measured is for reference only.
 5. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
 6. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.

7. Do NOT use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
8. Do NOT under any circumstances use the device on neonates or infants.
9. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
10. Do NOT use the equipment in places where aerosol sprays are being used or where oxygen is being administered.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

BEFORE YOU BEGIN	09
Important Information	09
Intended Use	10
Test Principle	10
Contents of System	11
Meter Overview	12
Display Screen	13
Test Strip	14
Setting the Meter	15
THE FOUR MEASURING MODES	18
Control Solution Testing	19
Performing a Control Solution Test	19
TESTING WITH BLOOD SAMPLE	22
Preparing the Puncture Site	23
Performing a Blood Glucose Test	24
METER MEMORY	27
Reviewing Test Results	27
Reviewing Blood Glucose Day Average Results	28
Downloading Results onto a Computer	29
MAINTENANCE	30
Battery	30
Caring for Your Meter	31
Disinfecting Procedures	32
Caring for Your Test Strips	34
Important Control Solution Information	35
SYSTEM TROUBLESHOOTING	36
Result Readings	36
Error Messages	37
Troubleshooting	39
DETAILED INFORMATION	40
SPECIFICATIONS	42

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values.
- Use only fresh whole blood samples to test the patient's blood glucose. Using other substances will lead to incorrect results.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult the healthcare professional before use.
- Critically ill patients should not be tested with a blood glucose monitoring system.

Intended Use

This system is intended for use outside the body (*in vitro* diagnostic use) by professionals for multiple-patient use in professional healthcare settings as an aid to monitoring the effectiveness of diabetes control. It is intended to be used for the quantitative measurement of glucose (sugar) in fresh whole blood samples. It is intended for use by professionals testing with capillary or venous and samples.

It should not be used for the diagnosis of or screening for diabetes, or testing on newborns.

The system is only used with single-use, auto-disabling lancing devices. It can only be used with FORA GD50 test strips and FORA control solutions.

Test Principle


Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

This system uses enzyme of Glucose Dehydrogenase (GDH) from *E. coli*.

Contents of System

The **FORA** GD50 system kit includes:

1

2

3

4

5

- | | |
|--------------------------------|----------------------------|
| 1 Meter | 4 Daily Log Book |
| 2 Owner's Manual/Warranty Card | 5 One AAA alkaline battery |
| 3 Protective Case | |


Test strips, control solutions and sterile lancets may be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

The USB cable is not provided in the kit but an accessory for GD50 which provides data transmission feature connecting to your PC (see section Downloading Results onto a Computer). Please contact local customer services or the place of purchase for obtaining the cable.

NOTE:

If any items are missing from your kit or opened prior to use, or to obtain test strips and control solutions, please contact the place of purchase or call the Customer Service for assistance.

Meter Overview

1 Display Screen

2 Test Strip Ejector

Eject the used strip by pushing this button.

3 Test Slot

Insert test strip here to turn the meter on for testing.

4 SET Button

Enter and confirm the meter settings.

5 Data Port

Download test results with a cable connection.


6 M Button

Display the meter memory; silence a reminder alarm; change a setting in SET mode.

7 Battery Compartment

8 Strip Indication Light

Display Screen

- | | |
|-----------------------------|---------------------------------|
| 1 Blood Drop Symbol | 9 Measuring Modes |
| 2 Test Strip Symbol | 10 Date |
| 3 Low Battery Symbol | 11 Control Solution Mode |
| 4 Test Result | 12 Memory Symbol |
| 5 Measurement Unit | 13 Error Message |
| 6 Reminder Alarm | 14 Ketone Warning |
| 7 Day Average | |
| 8 Time | |

Test Strip

Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

ATTENTION:

The front side of test strip should face up when inserting test strip.


Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE:

The **FORA GD50** monitor should only be used with **FORA GD50** Test Strips. Using other test strips with this meter can produce inaccurate results.

Setting the Meter

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

Entering the Setting Mode

Start with the meter off (no test strip inserted). **Press SET.**

1. Setting the date

With the year flashing, press **M** until the correct year appears. Press SET.

With the month flashing, press **M** until the correct month appears. Press SET.

With the day flashing, press **M** until the correct day appears. Press SET.

2. Setting the time format


Press and release **M** to select the desired time format — 12h or 24h. Press SET.

3. Setting the time

With the hour flashing, press **M** until the correct hour appears. Press SET.



With the minute flashing, press **M** until the correct minute appears. Press SET.

4. Deleting the memory

With “dEL” and a flashing “**M**” on the display, press **M** and select “no” to keep the results in memory then press SET to skip.

To delete all the results, press **M** and select “yes” and “**M**” are displayed on the meter, press SET to delete the memory.

5. Setting the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays “On” or “OFF” and . Press **M** to turn on or turn off to set the first reminder alarm.

Press **M** to select “On”, then press SET to set the hour. When the hour is flashing, press **M** to add an hour. Press SET to confirm and go to minutes, press

M to add one minute. Hold **M** longer to add faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.

If you want to turn off an alarm, find the alarm number by pressing SET in the setting mode, press **M** to change from “On” to “OFF”.

At the time of your alarm, the meter will beep and automatically turn on. You can press **M** to silence the alarm and insert a test strip to begin testing. If you do not press **M**, the meter will beep for 2 minutes then switch off. If you do not want to test at this time, press **M** to switch off the meter.

Congratulations! You have completed all settings!

NOTE:


- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

THE FOUR MEASURING MODES



The meter provides you with four modes for measuring, General, AC, PC and QC.

MODES	USE WHEN
General (displays as “Gen”)	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution



You can switch between each mode as follows:

1. Start with the meter switched off.
Insert a test strip to turn on the meter.
The screen will display a flashing “” and “Gen”.

2. Press **M** to switch between General, AC, PC and QC mode.

Control Solution Testing


Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

Do a control solution test when:


- you first receive the meter,
- at least once a week to routinely check the meter and test strips,
- you begin using a new vial of test strips,
- you suspect the meter or test strips are not working properly,
- patients' blood glucose test results are not consistent with how they feel, or if you think the results are not accurate,
- practicing the testing process, or you have dropped or think you may have damaged the meter.

Performing a Control Solution Test

1. Insert the test strip to turn on the meter

Insert the test strip into the meter. Wait for the meter to display “” and “”.



2. Press M to mark this test as a control solution test

With “**QC**” displayed, the meter will store your test result in memory under “**QC**”. If you press **M** again, the “**QC**” will disappear and this test is no longer a control solution test.


WARNING:

When doing the control solution test, you have to mark it so that the test result will **NOT** mix with the blood glucose **TEST RESULTS** stored in the memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

3. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

4. Read and compare the result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

(136 mg/dL = 7.6 mmol/L; 109-165 mg/dL = 6.1-9.2 mmol/L)

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Please contact the place of purchase or call the Customer Service for assistance.

NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for human blood glucose level.
- See the **Meter Maintenance** section for important information about your control solutions.

TESTING WITH BLOOD SAMPLE

WARNING: Potential Biohazard!

Healthcare professionals using this system on multiple patients should follow the infection control procedure approved by their facility. All products or objects which come in contact with human blood, even after cleaning, should be handled as if capable of transmitting viral diseases.

Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. Please follow the infection control procedure of lancing device approved by your facility.

To reduce the chance of infection:

- Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples.
- Wear a new pair of clean gloves before testing each patient. Change gloves between patients.
- Wear protective glasses and/or other protective clothing if necessary.
- Never share a lancet or the lancing device.
- Only auto-disabling, single use lancing devices may be used with this device.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

We recommend you perform disinfection procedures between each patient. Please refer to the section Cleaning and Disinfection Procedures for complete instructions. After disinfection, used gloves should be removed and hands washed before proceeding to the next patient.

Please refer to the following practice guidelines for more information about the correct procedure:

- Biosafety in Microbiological and Biomedical Laboratories (BMBL) found at <http://www.cdc.gov/biosafety/publications/bmbl5/>

- “Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Third Edition” Clinical and Laboratory Standards Institute (CLSI) M29-A3.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from an alternate site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:


- Wash and dry your hands before starting.
 - Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
 - Rub the puncture site for about 20 seconds before penetration.
 - Use a clear cap (included in the kit) while setting up the lancing device.
-
- **Fingertip testing**

Press the lancing device's tip firmly against the lower side of the patient's fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Performing a Blood Glucose Test




1. Insert the test strip to turn on the meter

Wait for the meter to display “ ” and “ ”.

2. Select the appropriate measuring mode by pressing M

For selecting the measurement mode, please refer to the “FOUR MEASURING MODES” section.

3. Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. Wipe off the first appeared drop of blood with a clean cotton swab. The size of the drop should be at least as big as ● (actual size), which is **0.5** microliter (μL) of volume. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

4. Apply the sample

Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

NOTE:


- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

5. Read your result

The result of the blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.

(92 mg/dL = 5.1 mmol/L)

6. Eject the used test strip

To eject the test strip, point the strip at a disposal container for sharp items. The meter will switch itself off automatically after the test strip is ejected.

WARNING:



The used lancet and test strip may be biohazardous . Please discard them carefully according to your local regulations.

Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination. For more information, please refer to the “Cleaning and Disinfection” section.

METER MEMORY

The meter stores the 450 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, **start with the meter switched off.**


Reviewing Test Results

1. Press and release M

“” will appear on the display. Press **M** again, and the first reading you see is the last blood glucose result along with date, time and the measurement mode.


(300 mg/dL = 16.7 mmol/L)



2. Press M to recall the test results stored in the meter each time you press.

After the last test results, press M again and the meter will be turned off.

(200 mg/dL = 11.1 mmol/L)

Reviewing Blood Glucose Day Average Results

1. Press and release **M**

When “**M**” appears on the display, keep pressing **M** for 3 seconds until the flashing “**DAY AVG**” appears. Release **M** and then your 7-day average result measured in general mode will appear on the display.

(114 mg/dL = 6.3 mmol/L)

2. Press **M to review** 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

3. Exit the meter memory. Keep pressing **M** and the meter will turn off after displaying the last test result.

NOTE:

- Any time you wish to exit the memory, hold **M** for 5 seconds or leave the meter idle for 3 minutes. The meter will turn off automatically.
- Control solution results are **NOT** included in the day average.
- If using the meter for the first time, “- - -” displays when you recall the test results or review the average result. It indicates that there is no test result in the memory.


Downloading Results onto a Computer

Data Transmission Via Cable

You can use the meter with an interface cable and the Health Care Software System to view test results on your personal computer. To learn more about the Health Care Software System or to obtain an interface cable separately, please contact local customer services or the place of purchase for assistance.

1. Obtaining the required cable and installing the software

To download the Health Care Software System, please visit the ForaCare website: www.foracare.com

2. Connecting to a personal computer

Connect the cable to a serial port on your computer. With the meter switched off, connect the other end of the interface cable to the meter data port. “PC” will appear on the meter display, indicating that the meter is in communication mode.

3. Data transmission

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

WARNING:

- While the meter is connecting to the PC, it will be unable to perform a blood glucose test.

MAINTENANCE

Battery

Your meter comes with one 1.5V AAA size alkaline battery.

Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low.


The “” appears with **E-b, Error and LOW:**

The power is not enough to do a test. Please change the batteries immediately.

Replacing the Battery

To replace the batteries, make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old battery and replace with one 1.5V AAA size alkaline battery.
3. Close the battery cover. If the battery is inserted correctly, you will hear a “beep” afterwards.

NOTE:

- Replacing the battery does not affect the test results stored in the memory.

- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Battery might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Caring for Your Meter

Wash and dry your hands thoroughly before use or handling of the meter and strips.

What is Cleaning and Disinfection?

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter?

Clean the meter when you see any dirt on it. You should disinfect the meter between each patient to prevent infection.

How to clean and disinfect the meter?

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt, blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. **Do NOT use organic solvents to clean the meter.**


We recommend using the product listed below to clean and disinfect your meter. It has been shown to be safe when used on the **FORA GD50** blood glucose meter.

► **Micro-Kill+™ (Micro-Kill Plus™) by Medline**
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

- 1.** Put on non-sterile gloves.
- 2.** Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

- 3.** Wipe all exterior surface of the meter including the display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot, be careful not to allow excess liquid to get inside. Keep the meter moist with disinfection solution contained in the wipe for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe. Use two or more wipes if necessary.

- 4.** Remove the wipe. Allow the meter surface to dry completely.
- 5.** Discard the used wipes and never reuse them.
- 6.** Remove and discard gloves in appropriate receptacles and wash hands.

Each cleaning and disinfection cycle includes a pre-cleaning step with one wipe and a disinfection step with a second wipe.

This device has been validated to withstand 10,000 cleaning and

disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 5 cleaning and disinfection cycles per day over 5 years, the expected life of device. The meter should be replaced after the validated number of disinfection cycles or the warranty period, whichever comes first.

Improper system cleaning and disinfection may result in meter malfunction. Stop using the meter if you see any signs of deterioration such as the LCD display is cracked or cloudy, the buttons no longer function, or the outer casing is cracked. Please contact customer service at 1-888-307-8188 for a replacement meter if any of the signs of deterioration are noticed.

For more information regarding the risk of transmission of bloodborne pathogens, please refer to:

- *“FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication” (2010)*

<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

- *“CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010)*

<http://www.cdc.gov/injectionsafety/Fingertick-DevicesBGM.html>

NOTE :

- Do **NOT** clean and disinfect the meter while performing tests.
- Please follow the instructions on the package label of Micro-Kill+™ disinfecting wipes for safe use of the wipes.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot, the battery compartment or the eject button.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.

- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

Meter Storage

- Storage conditions: -4 °F to 140 °F (-20 °C to 60 °C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage conditions: 35.6 °F to 89.6 °F (2 °C to 32 °C), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.

- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 68°F to 77°F (20°C to 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6°F to 86°F (2°C to 30°C). Do **NOT** freeze.

SYSTEM TROUBLESHOOTING


If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair yourself and never try to disassemble the meter under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
 The image shows a digital display with 'Lo' in large characters. Below it, in smaller text, is 'Gen', '8-10', and '10:00 AM'.	Appears when your blood glucose test result is below the lower measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
 The image shows a digital display with '300' in large characters. Below it, in smaller text, is 'KETONE? mg/dL', 'Gen', '8-10', and '10:00 AM'.	Appears when your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
 The image shows a digital display with 'Hi' in large characters. Below it, in smaller text, is 'Gen', '8-10', and '10:00 AM'.	Appears when your blood glucose test result is higher than the upper limit of measurement, which is 600 mg/dL (33.3 mmol/L).

Error Messages

MESSAGE	WHAT IT MEANS	WHAT TO DO
	Appears when the battery cannot provide enough power for a test.	Replace the battery immediately.
	Appears when a used test strip is inserted	Repeat with a new test strip.
	Appears when test strip is removed while counting down.	Re-test with a new test strip.
	Appears when ambient temperature is below the system operation range.	System operation range is 50°F to 104°F (10°C to 40°C). Repeat the test after the meter and test strip are in the above temperature range.
	Appears when ambient temperature is above the system operation range.	

MESSAGE	WHAT IT MEANS	WHAT TO DO
  Code  	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact your customer services for assistance.

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the battery.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer services.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 68 °F to 77 °F (20 °C to 25 °C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

DETAILED INFORMATION

Reference values

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining blood glucose levels close to normal can reduce the risk of diabetes complications by up to 60%*¹. The results provided by this system can help the healthcare professionals monitor and adjust treatment plan to gain better control of diabetes. The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2012). Clinical Practice Recommendations. Diabetes Care, 35 (Supplement 1): S1–S100.

*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Comparing Meter and Laboratory Results

The meter provides you with plasma blood equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results may be affected by factors and conditions that do not affect laboratory results in the same way. To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before going to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before doing comparison tests, if possible.

- Take your meter with you to the lab.

While staying at the lab:

Make sure that the samples for both tests are taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray-top test tube.
- Use fresh capillary or venous blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten, exercised, taken medication or experienced stress^{*2}. In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test^{*3}. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a meter result to be different from a laboratory result.

References

*2: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*3: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

SPECIFICATIONS

Model No.: FORA GD50

Dimension & Weight: 88 x 52 x 14 mm, 49 g

Power Source: one 1.5V AAA alkaline battery

Display: LCD

Memory: 450 measurement results with respective date and time

External Output: RS232 PC interface

Expected Service Life: 5 years

Auto electrode insertion detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature Warning

Operating Condition:

50 °F to 104 °F (10 °C to 40 °C), below 85% R.H. (non-condensing)

Storage/Transportation Conditions:

-4 °F to 140 °F (-20 °C to 60 °C), below 95% R.H.

Strip Storage/Transportation Conditions:

35.6 °F to 89.6 °F (2 °C to 32 °C), below 85% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

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