FORA® P30 Plus Series

Blood Pressure Monitoring System

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Manual del Propietario


TABLE OF CONTENTS

Important Safety Precautions	4
Before You Begin	4
Intended Use	4
Test Principle	4
Contents of System	5
Meter Overview	6
Display Screen	7
Setting the meter	7
Testing Your Blood Pressure	9
Fitting the Cuff Properly	10
Proper measurement position	11
Taking a Single Measurement	11
Averaging Measurement Mode	12
Auscultatory Mode	13
Recalling Past Readings	14
Reviewing Test Results	14
Data Transmission (for P30 Plus BLE only)	16
Data Transmission via Bluetooth	16
Maintenance	17
Battery	17
Using AC Adapter (Optional)	18
Caring for Your Meter	19
System Troubleshooting	19
Error Messages	19
Troubleshooting	20
Reference Values	20
Federal Communications Commission (FCC)	
Statement	22
Specifications	26

Important Safety Precautions

Read Before Use

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or damaged.
4. Do **NOT** use the equipment where aerosol sprays are being used, or where oxygen is being administered.
5. Do **NOT** use under any circumstances on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured are for reference only. Always consult your doctor to have the results interpreted.
7. Keep the equipment and its flexible cord away from hot surfaces.
8. Do **NOT** apply the cuff to areas other than the place directed.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Before You Begin

Intended Use

The **FORA P30 Plus Series** Blood Pressure Monitoring System (Model: P30 Plus / P30 Plus BLE) is a device intended to measure non-invasively the systolic and diastolic blood pressure and pulse rate of an adult individual. The blood pressure is measured by using an inflatable cuff wrapped around the arm. This system should be used for the testing on people over age of 18.

The device detects the appearance of irregular heartbeat during measurement, and gives a warning signal with the reading once the irregular heartbeat is detected.

Test Principle


Blood pressure is measured non-invasively at the arm based on oscillometric method.

This device is **NOT** able to take measurements in the presence of

common arrhythmia, such as atrial or ventricular premature beats or atrial fibrillation. It may produce reading error.

Contents of System

The **FORA P30 Plus Series** system kit includes:

1 Meter

2 Owner's Manual

3 Protective Wallet

4 Quick Start User Guide /
Warranty Card



5 Four AA Alkaline Batteries

6 Pressure Cuff

NOTE:

- If any items are missing from your kit or opened before use, please contact local customer service or place of purchase for assistance.

Meter Overview

Pressure Cuff

- 1. Display Screen**
- 2. M Button**
To enter the meter memory / To start the wireless data transmission (P30 Plus BLE only)
- 3. On/Off Button**
To start a single non-invasive blood pressure measurement.
- 4. AVG / S Button**
Blood pressure average mode / Enter and confirm the meter setting.
- 5. Bluetooth Indicator (P30 Plus BLE only)**
- 6. Battery Compartment**
- 7. AC Adapter Port**
Connect to a power supply.
- 8. Air Jack**
- 9. Pressure Cuff**
- 10. Air Tube**
- 11. Air Plug**
Connect to air jack.

Display Screen

- | | |
|---|-------------------------------|
| 1. Date | 7. Pulse Rate |
| 2. Systolic Pressure Value | 8. Day Average Results |
| 3. Diastolic Pressure Value | 9. Day / Night Time Symbol |
| 4. Measurement Unit | 10. Diastolic Pressure Symbol |
| 5. Memory Mode Symbol | 11. Battery Symbol |
| 6. IHB Detection Symbol
(Irregular Heart Beat) | 12. Systolic Pressure Symbol |
| | 13. Time |




Setting the meter

Before using your meter for the first time or if you change the battery, you should check and update these settings. Make sure you complete the steps below and save the desired settings.

Enter the Setting Mode



Start with the meter off. Press **AVG / S** until the meter turns on. The following settings could be set up in sequence:



1. Set the Date

With the YEAR / MONTH / DAY flashing in sequence, press **M** to select the number. Press **AVG / S** to confirm.

2. Set the Time & Time Format

With the HOUR / MINUTE flashing in sequence, press **M** to select the number. Press **AVG / S** to confirm.



With the HOUR and MINUTE flashing at the same time, press **M** to select the time format. Press **AVG / S** to confirm.


3. Delete the Memory

With “dEL” and “ ” on the display, keep all the results by pressing **AVG / S**. To delete all the results, press **M** twice.

4. Bluetooth Pairing (for P30 Plus BLE only)

With “PAR” on the display, press **M** to select “On”, and “PCL” on display with blinking Bluetooth indicator. Then complete the Bluetooth pairing (see the section “**Data Transmission via Bluetooth**”).

When the setting is complete, the meter will switch off automatically.

NOTE:

- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

Congratulations! You have completed all settings!

Testing Your Blood Pressure


This meter provides you different ways to measure your blood pressure. Select from the operations below:

- **Single measurement**
Perform an individual blood pressure measurement.
- **Averaging measurement**
Automatic perform three (3) consecutive blood pressure measurements and display the average result in the end.
- **Auscultatory measurement**
Only well trained persons may use this mode to measure blood pressure manually. The manual operation method involves applying a stethoscope and listening to the pulse during deflation while the air is slowly let out from the cuff (the Korotkoff method). Please ask your healthcare professional to train you how to use auscultatory mode.


Before Measurement

- Avoid caffeine, tea, alcohol and tobacco at least 30 minutes before measurement.
- Wait 30 minutes after exercising or bathing before measurement.
- Sit or lie down for at least 10 minutes. Take 5-10 minutes break between measurements. This break can be longer if necessary, depending on your physical conditions.
- Do not measure when feeling anxious or tense.
- Keep the records for your doctor as reference.
- Blood pressure varies between each arm. Always measure your blood pressure on the same arm.


Fitting the Cuff Properly

1. Connect the air plug of the tubing to the air jack of the meter.


2. Assemble the cuff as shown. The smooth surface should be inside the cuff loop and the metal D-ring should not touch the patient's skin.

3. Stretch the left (or right) arm in front with the palm facing up. Slide and place the cuff onto the arm with the air tube and artery mark region (in red) toward the lower arm.

Wrap and tighten the cuff above the elbow.

The red line on the edge of the cuff should be approximately **0.8 to 1.2 inches (2 to 3 cm)** above the elbow. Align the tube over the main arteries on the inside.

4. Leave a little free space between the arm and the cuff; you should be able to fit 2 fingers between them. Clothing must not restrict the arm. Remove all clothing covering or constricting the measurement arm.

5. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around the upper arm.

Proper measurement position

1. Sit down for at least 10 minutes before measuring.
2. Place your elbow on a flat surface. Relax your hand with the palm facing up.
3. **Make sure the cuff is about the same height as the location of your heart.**
Press **On/Off** button. Remain still and do not talk or move during the measurement.
4. Measurement is in progress.

After the meter is turned on, the cuff will begin to inflate automatically.

NOTE:

- If the cuff is relatively lower (higher) than the heart, the obtained blood pressure value could be higher (lower) than the actual value. A 15 cm difference in height may result in an error around 10 mmHg.
- Blood pressure readings may be affected if you do not use in proper measurement position.

Taking a Single Measurement

Always apply the pressure cuff before turning on the meter.

1. Press **On/Off** button to turn the meter on. Then the cuff will begin to inflate automatically. The heart symbol "♥" will flash when a pulse is detected during the inflation.




2. After the measurement, the meter displays the systolic pressure, diastolic pressure and pulse rate.
3. Press **On/Off** button to turn the meter off.

NOTE:

- If you press **On/Off** button during measurement, the meter will be turned off.
- If the pulse rate symbol is shown as “” instead of “”, this indicates that the meter has detected an irregular heartbeat (IHB). Please consult your doctor if you have question about IHB.

Averaging Measurement Mode

Always apply the pressure cuff before turning on the meter.

1. Press **AVG / S** button to turn the meter on. Then the cuff will begin to inflate automatically.
2. After the first measurement is complete, the meter will start counting down before the second measurement begins. The number on the right represents the remaining countdown between each measurement. The device will take three (3) measurements consecutively with an interval of 20 seconds.

Note: When the difference between the first and second systolic pressure is more than 15mmHg, the time interval for third measurement will be 40 seconds.


3. After taken three measurements, the results are averaged to produce the blood pressure measurement with “AVG” symbol display on the screen.
4. Press **On/Off** button to turn the meter off.

Auscultatory Mode

Only well trained persons may use this mode to measure blood pressure manually. This manual method involves applying a stethoscope to the arm and listening to the pulse while the air is slowly let out from the cuff (the Korotkoff method). Please ask your health-care professional to train you how to use auscultatory mode.

The systolic pressure is the maximum pressure in an artery at the moment when the heart is beating and pumping blood through the body. The diastolic pressure is the lowest pressure in an artery in the moments between beats when the heart is resting.

1. Place a stethoscope on the patient's arm where there is a pulse. Apply the cuff to the patient's upper arm and hold in place with Velcro.
2. Press and hold **On/Off** button until "AUS" on the display. An AU symbol will appear on the display with a beep. Release the button. Then cuff begins to inflate automatically. After reaching the cuff pressure, the deflation begins.

3. You can adjust the inflation pressure by keep pressing **M** to inflate if necessary. If you release **M**, the inflation will stop and deflation begins.
4. The systolic pressure is measured when the operator first hears the pulse. Take the reading according to the value on the display. This reading indicates the systolic pressure. You have to mark the reading due to the result will not be stored in the memory.
5. The diastolic pressure is measured from the moment the operator is unable to hear the sound of the pulse. Take the reading according to the value on the display. This reading indicates the patient's diastolic pressure. You have to mark the reading due to the result will not be stored in the memory.

6. The meter will return to the screen of 0 mmHg after the measurement is completed. To turn off the meter, press and hold **On/Off** button. Or it will switch off automatically after being idle for 3 minutes.

Recalling Past Readings

Your meter stores the 200 most recent test results with date and time in its memory. To enter the meter memory, start with the meter off.


Reviewing Test Results

For P30 Plus:


Reviewing Single Measurement Results

1. Press **M** to enter memory mode.

2. With "01" and "  " on the display, and the first reading you see is the last blood pressure result along with date and time.

3. Press **M** to review all test results stored in the meter.
4. Press **On/Off** button to turn the meter off.

Reviewing Day Average Results

1. Press and hold **M** to enter day average memory mode.
With "DAY AVERAGE" and "  " on the display, and the first reading you see is the day average result measured in general mode will

appear on the display.

2. Press **M** to review the Day-time (☀️) average result.
3. Press **M** to review the Night-time (🌙) average result.
4. Press **On/Off** button to turn the meter off.

For P30 Plus BLE:


Reviewing Day Average Results (1-3)


1. Press **M** to enter memory mode. The day average result measured in general mode will appear on the display.
2. Press **M** to review the Day-time (☀️) average result.
3. Press **M** to review the Night-time (🌙) average result.

Reviewing Single Measurement Results (4-5)

4. Press **M** to scroll through all readings stored in the meter.
5. Press **On/Off** button to turn the meter off.

NOTE:

- If there are no records in memory, “---” displays when you recall the

test results or review the average results. This indicates that there is no test result in the memory.

- ☼ Day-time average means the average records for measurements taken between 4:00 A.M. and 11:59 A.M.
- 🌙 Night-time average means the average records for measurements taken between 6:00 P.M. and 11:59 P.M.

Data Transmission (for P30 Plus BLE only)

FORA P30 Plus Series has 2 versions: P30 Plus (no data transmission function) and P30 Plus BLE (Bluetooth Connectivity). Please check your meter box for your meter version.

Data Transmission via Bluetooth

You can transmit your data from the meter to your devices (e.g. smart phone, tablet, PC...) via Bluetooth. Please contact your local customer service or place of purchase for assistance. Please note that you must pair your meter to your device before transmitting data.

Pairing with your mobile device

1. Press and hold **M** until "PCL" on the display with flashing Bluetooth indicator.
2. On your mobile device, turn on the Bluetooth function and follow the instruction on your APP to pair the device (e.g., scan for the meter to add it into APP).
3. After successfully pairing the APP to your meter, the Bluetooth function on the meter will then start to transmit your data to your APP.

Bluetooth indicator on the meter

BLUETOOTH INDICATOR	STATUS
Flashing Blue	The Bluetooth function is enabled and waiting for connection.
Solid Blue	The Bluetooth connection is established.

NOTE:

- While the meter is in transmission mode, it will be unable to perform a test.

- The meter is compatible with the following devices, iOS (6 or above) and Android system (4.3 or above). Make sure device is in proper receiving range of the meter before transmitting data.
- A compatibility issue between your mobile device and the meter might occur due to the way Bluetooth functionality is implemented by the various mobile device manufacturers.


Maintenance

Battery

Your meter comes with four 1.5V AA size alkaline batteries.

Low Battery Signal

The meter will display the two messages below to alert you when the meter power is getting low.

1. The “” appears along with a displayed message: The meter is functional and the result remains accurate, but it is time to change the battery.
2. The “” appears with E-b: There is not enough power to do a test. You must change the batteries immediately.

Replacing the Battery

To replace the battery, make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old batteries and replace with four 1.5V AA size alkaline batteries.

3. Close the battery cover.

NOTE:


- Replacing the batteries does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Using AC Adapter (Optional)

Connect an AC adapter to the meter.


1. Connect AC adapter plug to AC adapter jack of the meter.

2. Plug AC adapter power plug into an electrical outlet. Press **On/Off** button to start the measurement.

Remove an AC adapter from the meter.

1. **When the meter is off**, remove the AC adapter power plug from the electrical outlet.
2. Disconnect the AC adapter plug from the AC adapter jack of the meter.

Caring for Your Meter

To avoid the meter and pressure cuff attracting dirt, dust or other contaminants, please wash hands thoroughly with soap and water before use.

Cleaning (after use or transfer the device to another patient's use)

- To clean the meter exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do **NOT** flush with water.
- Do **NOT** use organic solvents to clean the meter.
- Do **NOT** wash the pressure cuff.
- Do **NOT** iron the pressure cuff.

Meter Storage

- Storage condition: -13°F to 158°F (-25°C to 70°C), 10% to 93% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

System Troubleshooting

If you follow the recommended action but the problem persists, or error messages occur, please call your local customer service. Do not attempt to repair yourself and never try to disassemble the meter under any circumstances.

Error Messages

MESSAGE	WHAT IT MEANS	WHAT TO DO
	Big arm movement or vibration and unable to detect the pulse during the measurement.	The meter will have a second measurement automatically. Try not to move the arm and remain still during measurement.
E-1 E-2	Inflation / deflation or pressure error.	Refit cuff tightly and correctly. Relax and repeat the measurement. If the error still remains, please contact local customer service for help.
E-4 E-5 E-6	Blood pressure measurement error.	

E-A E-E	Problems with the meter.	Please contact local customer service for help.
E-b	Battery is too low.	Replace the batteries or input AC adapter.

Troubleshooting

1. If nothing is displayed after pressing On/Off button:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

2. If the heart rate is higher/lower than user's average:

POSSIBLE CAUSE	WHAT TO DO
Movement during measurement.	Repeat measurement.
Measurement taken just after exercise.	Rest at least 30 minutes before repeating measurement.

3. If blood pressure is higher/lower than user's average measurement:

POSSIBLE CAUSE	WHAT TO DO
May not be in correct position while measuring.	Adjust to the correct position to measure.
Blood pressure naturally varies from time to time.	Keep in mind for the next measurement.

4. If the cuff inflates again while measuring

POSSIBLE CAUSE	WHAT TO DO
Cuff is not fastened.	Fasten the cuff again.
If user's blood pressure is higher than the pressure the device has inflated, the device will automatically increase the pressure and start to inflate again. Stay relaxed and wait for the measurement.	

Reference Values

Clinical studies show that adult diabetes is often accompanied by

elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment^{*1}. Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

The recommended blood pressure range is as below:

Classification	Systolic Pressure (mmHg)	Diastolic Pressure (mmHg)
Hypotension ^{*2}	< 90	< 60
Normal ^{*3}	< 120	< 80
Elevated ^{*3}	120 – 129	< 80
Stage 1 Hypertension ^{*3}	130 – 139	80 – 89
Stage 2 Hypertension ^{*3}	≥ 140	≥ 90
Hypertensive Crisis ^{*3}	> 180	> 120

*1: National Heart, Lung, and Blood Institute, Diseases and Conditions.

*2: The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure. NIH Publication. 2003. No. 03-5233.

*3: New ACC (American College of Cardiology) and AHA (American Heart Association) guidelines for the detection, prevention, management and treatment of high blood pressure (2017).

Federal Communications Commission (FCC) Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference and
2. This device must accept any interference received, including interference that may cause undesired operation of the device.

15.105(b)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

FCC RF Radiation Exposure Statement:

1. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. For portable operation, this device has been tested and meets FCC RF exposure guidelines. When used with an accessory that contains metal may not ensure compliance with FCC RF exposure guidelines.

Manufacturer's declaration-electromagnetic emissions		
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance (for home and professional healthcare environment)
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Manufacturer's declaration-electromagnetic immunity			
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home and professional healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient /burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0% UT; 0.5 cycle 0% UT; 1 cycle 70% UT; 25/30 cycles Voltage interruptions: 0% UT; 250/300 cycle	Voltage dips: Not applicable Not applicable Not applicable Voltage interruptions: Not applicable	Mains power quality should be that of a typical home and professional healthcare environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical home and professional healthcare environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Manufacturer's declaration-electromagnetic immunity

The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below.
The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home and professional healthcare environment)
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	Not applicable Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	10 V/m: 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m: 80 MHz – 2,7 GHz 80% AM at 1 kHz	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Recommended separation distance between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment (for home and professional healthcare) in which radiated RF disturbances are controlled.
The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Manufacturer's declaration-electromagnetic immunity							
Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment							
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.							
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	Compliance LEVEL (V/m) (for home and professional healthcare)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine	2	0,3	28	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9	9
745							
780							
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28	28
870							
930							
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28	28
1845							
1970							
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9	9
5500							
5785							
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.							
^{a)} For some services, only the uplink frequencies are included.							
^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal.							
^{c)} As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.							

Manufacturer's declaration-electromagnetic immunity Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields				
The device is intended for use in the electromagnetic environment (for home and professional healthcare) specified below. The customer or the user of the device should assure that it is used in such an environment.				
Frequencies	Test level [A/m]	Modulation	Dwell time [s]	Compliance LEVEL [A/m] (for home and professional healthcare)
30 kHz (a)	8	CW	3	8
134,2 kHz	65	Pulse modulation (b) 2,1 kHz	3	65 (c)
13,56 MHz	7,5	Pulse modulation (b) 2,1 kHz	3	7,5 (c)
Note:				
(a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME AND PROFESSIONAL HEALTHCARE ENVIRONMENT.				
(b) The carrier shall be modulated using a 50% duty cycle square wave signal.				
(c) r.m.s, before modulation is applied.				

Specifications

● System Performance

Power Source: Four 1.5V AA alkaline batteries or (optional) DC 6V power adapter (Input: AC 100 to 240V; Output: DC 6V, 1A, 6W)

Size of Meter w/o Cuff: 150 (L) x 100 (W) x 80 (H) mm

Weight: 300g approx. (without batteries)

Memory: Maximum 200 memory records

Power Saving: Automatic power off if system is idle for 3 minutes

Communication: Bluetooth (for P30 Plus BLE only)

System Operating Condition:

41°F to 104°F (5°C to 40°C), 15% - 93% R.H.

Meter Storage Condition:

-13°F to 158°F (-25°C to 70°C), 15% - 93% R.H.

Expected service life: 3 years (Meter) / 2 years (Cuff)

● Blood Pressure Measurement Performance

Measurement Unit: mmHg

Systolic Measurement Range: 60 to 255 mmHg

Diastolic Measurement Range: 30 to 195 mmHg

Pulse Rate Measurement Range: 40 to 199 beats / minute

Accuracy of Pressure: ± 3 mmHg or $\pm 2\%$ of reading

Accuracy of Pulse Rate: $\pm 4\%$ of reading

This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2.

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Distributed by Foracare, Inc. / Distribuido por Foracare, Inc.