
FORA[®]

Premium V12

Premium V12 Voice

Blood Glucose Monitoring System

Sistema de vigilancia de glucose sanguinea

Owner's Manual
Manual del Propietario

Dear FORA Premium V12 Series System Owner:

Thank you for purchasing the FORA Premium V12 Series Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read the following contents thoroughly and carefully.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the FORA Premium V12 Series Blood Glucose Monitoring System to easily monitor your blood glucose levels by yourself anywhere, any time.

If you have questions, just call ForCare Customer Care Service Center toll-free at 1-888-307-8188 (8:30 am - 5:00 pm PST, Mon. - Fri.). For assistance outside of these hours, please contact your healthcare professional. You can also visit www.foracare.com.

* Talking function is available on the FORA Premium V12 Voice System only. The voice output of FORA Premium V12 Voice System is available in English and Spanish.

Premium V12 Voice Model Only

IMPORTANT SAFETY PRECAUTIONS

READ BEFORE USE

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)
<http://wayback.archive-it.org/7993/20170111013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

1. Use this device **ONLY** for the intended use described.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only.
5. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
6. Use of this device in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
7. Do **NOT** use this device in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

TABLE OF CONTENTS

BEFORE YOU BEGIN	06
Important Information	06
Intended Use	07
Test Principle	07
Contents of System	08
Meter Overview	09
Display Screen	10
Test Strip	11
Setting the Meter	12
THE FOUR MEASURING MODES	16
Control Solution Testing	17
Performing a Control Solution Test	18
TESTING WITH A BLOOD SAMPLE	21
Preparing the Puncture Site	21
Alternative Site Testing	23
Performing a Blood Glucose Test	24
METER MEMORY	27
Reviewing Test Results	27
Reviewing Blood Glucose Day Average Results	28
Downloading Results on to a Computer	29
MAINTENANCE	30
Battery	30
Caring for Your Meter	31
Caring for Your Test Strips	34
Important Control Solution Information	35
SYSTEM TROUBLESHOOTING	36
Result Readings	36
Error Messages	37
Troubleshooting	39
DETAILED INFORMATION	40
SPECIFICATIONS	41

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than normal, and you do not have any symptoms of illness, repeat the test. If you have symptoms or continue to get results which are higher or lower than normal, consult with your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings lower than actual values may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult a healthcare professional before use.
- Critically ill patients should not be tested with a blood glucose monitoring system.

Intended Use

The **FORA Premium V12 Series** Blood Glucose Monitoring System is intended for use in the quantitative measurement of glucose in fresh capillary whole blood from the finger and the following alternative sites: palm, forearm and upper-arm. It is intended for use by people with diabetes mellitus at home as an aid in monitoring the effectiveness of diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and is not intended for use on neonates. It is intended to be used by a single person and should not be shared.

The alternative site testing in the **FORA Premium V12 Series** Blood Glucose Monitoring System can be used only during steady-state blood glucose conditions.


FORA Premium V12 Voice has some speaking functions but is not intended for use by the visually impaired.

Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Contents of System

Your new **FORA Premium V12** system kit includes:


- | | |
|----------------------------|---|
| 1 Meter | 4 Quick Start User Guide /
Warranty Card |
| 2 Owner's Manual | |
| 3 Protective Wallet | 5 Daily Log Book |
| | 6 2 x 1.5V AAA Batteries |

Test strips, control solution, sterile lancets and lancing device may be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

NOTE:


If any items are missing from your kit or opened prior to use, or to obtain test strips and control solutions, please contact the place of purchase or call Customer Service for assistance.

Meter Overview

- 1 Display Screen**
- 2 Test Strip Slot**
Insert test strip here to turn the meter on for testing.
- 3 SET Button**
Enter and confirm the meter settings.
- 4 Data Port**
Download test results with a cable connection.
- 5 M Button**
Enter the meter memory and silence a reminder alarm.
- 6 Strip indication light**

Display Screen

- | | |
|---|------------------------------|
| 1 Volume Symbol (For Premium V12 Voice Only) | 8 Time |
| 2 Test Result | 9 Day Average |
| 3 Ketone Warning | 10 Measurement Unit |
| 4 Memory Symbol | 11 Low Battery Symbol |
| 5 Control solution Mode | 12 Test strip symbol |
| 6 Date | 13 Blood Drop Symbol |
| 7 Test mode | |

Test Strip

Absorbent Hole

Touch to a drop of blood here, blood will be automatically drawn in.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Confirmation Window

Confirm if enough blood has been applied to the absorbent hole of the strip.

Contact Bars

Insert this end of the test strip into the test strip slot.

ATTENTION:

The front side of test strip should face up when inserting test strip.

Test results might be wrong if the contact bar is not fully inserted into the test slot.


NOTE:

The **FORA Premium V12 Series** Blood Glucose Monitoring System should only be used with **FORA Premium V12 Series** Test Strips. Using other test strips with this meter can produce inaccurate results.

Setting the Meter

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.


Entering the Setting Mode


Start with the meter off (no test strip inserted). **Press SET.**

1. Setting the date

With the year flashing, press M Button until the correct year appears. Press SET Button.

With the month flashing, press M Button until the correct month appears. Press SET Button.


With the day flashing, press M Button until the correct day appears. Press SET Button.

2. Setting the time

With the hour flashing, press M Button until the correct hour appears. Press SET Button.


With the minute flashing, press M Button until the correct minute appears. Press SET Button.

3. Deleting the memory

With “yes/no” flashing and “dEL” and “” symbol on the display, press M Button to switch decision.


If you would not delete the memory, press SET Button when flashing “no” on display.

When flashing “yes” on display, press SET Button to delete all internal records.

LCD will display “---” and “OK” after all records are deleted, press SET Button to continue the settings.


Premium V12 Voice model enters STEP 4 and STEP 5 for speaking volume setting and language setting mode; for Premium V12 model please go to STEP 6 directly.

4. Setting the speaking volume (For Premium V12 Voice only)


There are seven (7) speaking volume options to choose from. Press M Button until the desired speaking volume appears. To confirm your selection, press SET Button.

Volume 0 indicates that the speaking function is turned off, and “” will not display during testing. Volume 1 to 7 indicates speaking volume from low to high, and “” will be displayed throughout the testing.

5. Choosing a Language (For Premium V12 Voice only)

Press the M Button to select L1/L2. The default language for the meter is L1, which is English. To confirm your selection, press SET Button.


6. Setting the Reminder Alarm

You may set up any or all of the reminder alarms(1-4). The meter displays “On” or “OFF” and , press M Button to switch “On” or “OFF”. Press SET Button to set the first reminder alarm.

Press M Button to increase hour, then press SET Button to set the hour. Press M Button to increase minute, then press SET Button to set the minute and confirm. Go to the next alarm setting.

If you want to turn off an alarm, find the alarm number by pressing SET Button in the setting mode, press M Button to change from “On” to “OFF”.

At the time of your alarm, the meter will beep and automatically turn on. You can press M Button to silence the alarm or insert a test strip to begin testing.

If you do not press M Button, the meter will beep for 2 minutes then switch off automatically. If you do not want to test at this time, press M Button to switch off the meter.

Congratulations! You have completed all settings!

NOTE:


- These parameters can ONLY be changed in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

THE FOUR MEASURING MODES


The meter provides you with four modes for measuring, General, AC, PC and QC

Modes	Use when
General tests	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between each mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a flashing blood drop.

2. Press M Button to switch between General, AC, PC and QC mode

Control Solution Testing


Our control solutions contain a known amount of glucose that reacts with test strips and is used to ensure that the meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have all items prior to testing your blood glucose.

Test with a control solution when:

- you first receive the meter.
- you do your weekly routine quality testing of the meter and test strips.
- you begin using a new vial of test strips.
- you suspect the meter and or test strips are not working correctly.
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate.
- you are practicing the testing procedure.
- you have dropped or think you may have damaged the meter.

Performing a Control Solution Test


1. Insert the test strip to turn on the meter

Wait for the meter to display “” and “”.

**Premium V12 Voice model only.
Voice Guide:**

" Please apply blood into the strip."



2. Press M Button to mark this test as a control solution test

With “QC” displayed, the meter will not store your test result in memory. If you press the M Button again, the “QC” will disappear and this test is no longer a control solution test.

WARNING:

When doing the control solution test, you have to mark it so that the test result will not be stored in the memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

3. Apply control solution


Shake the control solution vial thoroughly before use.

Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap.

Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down.

To avoid contaminating the control solution, do not directly apply control solution onto a strip.

4. Read and Compare the Result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

Premium V12 Voice model only.
Voice Guide:

" The blood glucose is (number) milligrams per deciliter. "

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do NOT test your blood. Contact the local customer service or place of purchase for help.

NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **MAINTENANCE** section for important information about your control solutions.

TESTING WITH A BLOOD SAMPLE

Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets or the lancing device.
- Wash and dry your hands thoroughly before and after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “**MAINTENANCE**” section.
- If the meter is being operated by a person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by that person.

Sharing the lancing device and lancets may increase the risk of spreading infectious diseases. Lancing device must not be used on more than one person.

Preparing the Puncture Site


Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from an alternate site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.

Please follow the suggestions below before obtaining a drop of blood:


Wash and dry your hands before starting.

- Select the puncture site either at fingertips or another body part (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Rub the puncture site for about 20 seconds before penetration.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Use a clear cap while setting up the lancing device.

● **Fingertip testing**

Press the tip of the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger, a click indicates that the puncture is complete.

● **Blood from sites other than the fingertip**

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm or hand, avoid lancing the areas with obvious veins to avert excessive bleeding.

NOTE:


- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing

Important:

There are limitations with **Alternative Site Testing (AST)**. Please consult your health care professional before you perform **AST**.

- **Alternative site testing (AST) should not be used for calibrating continuous glucose monitors (CGMs).**
- **Alternative site testing (AST) should not be used for insulin dose calculations.**

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, **take a blood sample from your finger only.**


We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not match the way you feel.
- Your routine glucose results often fluctuate.

Performing a blood glucose test

1. Insert the test strip to turn on the meter

Wait for the meter to display the test strip  and flashing blood drop  symbols.



Premium V12 Voice model only.
Voice Guide:

" Please apply blood into the strip. "

2. Select the appropriate measuring mode by pressing M Button.

For selecting the measurement mode, please refer to the "THE FOUR MEASURING MODES" section.


3. Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. After penetration, discard the first drop of blood with a clean cotton swab. Gently squeeze the punctured area to obtain another drop of blood. Be careful NOT to smear the blood sample.


The volume of blood sample must be at least 0.7 microliter (μL) of volume.

( actual size).

4. Apply the sample

Hold the blood drop to touch the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the meter begins counting down. Do not remove your finger until the meter beeps.

NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

5. Read Your Result

The result of your blood glucose test will appear after the meter counts down to 0. This blood glucose result will automatically be stored in the memory.

**Premium V12 Voice model only.
Voice Guide:**

" The blood glucose is (number) milligrams per deciliter. "

6. Remove the used test strip and the lancet

Remove the test strip to a disposal container for sharp items. The meter will switch itself off automatically after the test strip is removed.

Always follow the instructions in the lancing device insert when removing the lancet.



WARNING:

The used lancet and test strip may be biohazardous. Please discard them properly.

METER MEMORY

The meter stores the 450 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, start with the meter switched off.

Reviewing Test Results


1. Press and release M Button, flashing “**M**” will appear on the display. Press M Button again, and the first reading you see is the last blood glucose result along with date, time and the measurement mode.

Premium V12 Voice model only.
Voice Guide:

"The blood glucose is (number) milligrams per deciliter."


2. Press the M button to recall the test results stored in the meter each time you press. After the last test result, press M again and the meter will be turned off.

Premium V12 Voice model only.
Voice Guide:


"The blood glucose is (number) milligrams per deciliter."

Reviewing Blood Glucose Day Average Results

1. Press and release M Button.

When “M” appears on the display, keep pressing M for 3 seconds until the flashing “DAY AVG” appears. Release M and then your 7-day average result measured in general mode will appear on the display.

2. Press M to review 14-, 21-, 28-, 60- and 90- day average results stored in memory.

3. Exit the meter memory. Press M after the 90-day average is displayed and the meter will turn off.

NOTE:

- Any time you wish to exit the memory, hold press **M** down for 5 seconds or leave the meter idle for 3 minutes. The meter will turn off automatically.
- Control solution results are **NOT** included in the day average.
- If there are no records in memory, “—” displays when you recall the test results or review the average result.


Downloading Results on to a Computer

Data Transmission Via RS232 Cable

You can use the meter with an RS232 cable and the Health Care Software System to view test results on your personal computer. To learn more about the Health Care Software System or to obtain an RS232 cable separately, please contact local customer services or the place of purchase for assistance.

1. Obtaining the required cable and installing the software

To download the Health Care Software System, please visit the ForaCare, Inc. website at www.foracare.com.

2. Connecting to a personal computer

Connect the cable to a RS232 port on your computer. With the meter switched off, connect the other end of the RS232 cable to the meter data port. "PC" will appear on the meter display, indicating that the meter is in communication mode.

3. Data transmission

To transmit data, follow the instructions provided with the software. Results will be transmitted with date and time. Remove the cable and the meter will automatically switch off.

WARNING:

While the meter is connecting to the PC, it will be unable to perform a blood glucose test.


MAINTENANCE

Battery

Your meter comes with two 1.5V AAA size alkaline batteries.

Low Battery Signal

The meter will display the message below to alert you when the meter power is getting low.


The  symbol appears with **E-b, Error and low:**

The power is not enough to do a test. Please change the batteries immediately.

Replacing the Battery

To replace the batteries, make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Align the positive and negative terminals of the batteries with the positive (+) and negative (-) signs engraved inside the battery compartment.
4. If the batteries are inserted correctly, you will hear a "beep" afterwards. Then close the battery cover.

NOTE:

- Replacing the batteries does not affect the test results stored in the memory.
- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Caring for Your Meter

Wash and dry your hands thoroughly before use or handling of the meter and strips.

What is Cleaning and Disinfection?

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter?

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter?

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. **Do NOT use organic solvents to clean the meter.**

We recommend using the product listed below to clean and disinfect your meter.

It has been shown to be safe when used on the FORA Premium V12 Series blood glucose meter.


► Micro-Kill+™ (Micro-Kill Plus™) by Medline
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact

Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all exterior surface of the meter including the display and buttons.

Hold the meter with the test strip slot pointing down and wipe the area around the test slot, be careful not to allow excess liquid to get inside. Keep the meter moist with disinfection solution contained in the wipe for a minimum of 2 minutes. Follow the instructions on the package label of disinfecting wipe. Use two or more wipes if necessary.

3. Remove the wipe. Allow the meter surface to dry completely.

4. Discard the used wipes.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-888-307-8188 for assistance.

This device has been validated to withstand 5,000 disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 1,000 tests per year on each meter and around 3 disinfection cycles per meter per day for 5 years, the expected life of device.

Stop using the meter if you see any signs of deterioration such as the LCD display is cracked or cloudy, the buttons no longer function, or the outer casing is cracked. Please contact the customer service at 1-888-307-8188 for a replacement meter if any of the signs of deterioration are noticed.

The meter should be replaced after the validated number of disinfection cycles or the warranty period, whichever comes first.

NOTE:

- Do **NOT** clean or disinfect the meter while performing a test.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to handling by an assistant.
- Do **NOT** allow cleaning or disinfecting solution to get in the test slot, battery compartment, or strip-ejection button.
- If you do get moisture in the test strip slot, wipe it away with the corner of a tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case and damage the electronic components or circuitry.

Meter Storage

- Storage conditions: -4°F to 140°F (-20°C to 60°C), 10% to 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage conditions: 39.2°F to 104°F (4°C to 40°C), 10% to 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months of first opening.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.


Important Control Solution Information

- Use only our control solutions with your meter.
- Do not use the control solution beyond the expiry date or 3 months after first opening. Write the opening date on the control solution bottle and discard the remaining solution after 3 months of first opening.
- It is recommended that the control solution test be done at room temperature (68°F to 77°F / 20°C to 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the bottle before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6°F and 86°F (2°C and 30°C). Do **NOT** freeze.


SYSTEM TROUBLESHOOTING


If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair the meter yourself and never try to disassemble the meter under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
 The image shows a digital display with 'Lo' in large digits. Above it is a small icon of a person with a plus sign. Below 'Lo' is the text 'Gen'. At the bottom, it shows '8-10 10:00 AM'.	Appears when your result is below measurement limit, which is less than 20 mg/dL.
The above symbol indicates hypoglycemia (low blood glucose). You should seek medical assistance immediately.	
 The image shows a digital display with '58' in large digits. Above it is a small icon of a person with a plus sign. Below '58' is 'mg/dL'. Below that is a small icon of a meter and the text 'Gen'. At the bottom, it shows '8-10 10:00 AM'.	Your actual reading appears when your glucose range from 20-240mg/dL
 The image shows a digital display with '300' in large digits. Above it is a small icon of a person with a plus sign. Below '300' is 'KETONE: mg/dL'. Below that is 'Gen'. At the bottom, it shows '8-10 10:00 AM'.	Ketone appears below the reading when your result is equal to or higher than 240 mg/dL and under 600 mg/dL. This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
 The image shows a digital display with 'Hi' in large digits. Above it is a small icon of a person with a plus sign. Below 'Hi' is the text 'Gen'. At the bottom, it shows '8-10 10:00 AM'.	Appears when your result is higher than the limit of measurement, which is higher than 600 mg/dL.

Error Messages

MESSAGE	WHAT SOUND DO YOU HEAR?	WHAT IT MEANS	WHAT TO DO
	Several beeps	Problem with the meter.	Review the instructions and repeat the test. If the meter still does not work, contact the local distributor for assistance.
	None		
	Several beeps		
	Several beeps	Appears when the batteries cannot provide enough power for a test.	Replace the batteries immediately.
	Several beeps	Problem with the meter.	Review the instructions and repeat the test. If the meter still does not work, contact the local distributor for assistance.
	Several beeps	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.

MESSAGE	WHAT SOUND DO YOU HEAR?	WHAT IT MEANS	WHAT TO DO
	Several beeps	Appears when test strip is removed while counting down, or other problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, contact the local distributor for assistance.
	Several beeps	Appears when a used test strip is inserted	Repeat with a new test strip.
	Several beeps	Appears when ambient temperature is below system operation range.	System operation range is 50° F to 104° F (10° C to 40° C). Repeat the test after the meter and test strip are in the operating temperature range.
		Appears when ambient temperature is above system operation range.	

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing  appears on the display.
Defective meter.	Please contact customer services.

3. If the control solution testing result is out of range.

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (68°F-77°F/20°C-25°C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

DETAILED INFORMATION

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association. Standards of Medical Care in Diabetes- 2018 Jan; 41(Supplement 1): S1-S2.

Please consult your doctor to determine a target range that works best for you.

SPECIFICATIONS

Model No.: Premium V12 / Premium V12 Voice

Dimension & Weight: 87.6 (L) x 54.5 (W) x 22 (H) mm, 47.8 g

Power Source: Two 1.5V AAA alkaline batteries

Display: LCD

Memory: 450 measurement results with respective date and time

External output: RS232

Expected service life: 5 years

Auto electrode insertion detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature Warning

Operating Condition: 50°F to 104°F (10°C to 40°C), 10% to 85% R.H.
(non-condensing)

Storage/Transportation Conditions: -4°F to 140°F (-20°C to 60°C),
10% to 95% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1,
IEC/EN 61326-2-6.

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1. "FDA Public Health Notification: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)

<http://wayback.archive-it.org/7993/20170111013014/http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

2. "CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010)

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Distributed by Foracare, Inc. / Distribuido por Foracare, Inc.