
FORA *TN'G*[®] Voice

Blood Glucose Monitoring System

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Owner's Manual
Manual del Propietario

Dear FORA TN'G Voice System Owner:

Thank you for purchasing the **FORA TN'G Voice** Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read and follow the instructions in this Owner's Manual.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the **FORA TN'G Voice** Blood Glucose Monitoring System to easily monitor your blood glucose levels.

If you have questions, just call Foracare Customer Care Service Center toll-free at 1-888-307-8188 (6:30am-5:30pm PST, Mon.-Fri.). You can also visit www.foracare.com.

IMPORTANT SAFETY PRECAUTIONS

READ BEFORE USE

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. Letter to Manufacturers of Blood Glucose Monitoring Systems Listed With the FDA

<https://www.fda.gov/medical-devices/industry-medical-devices/letters-industry>

2. Considerations for Blood Glucose Monitoring and Insulin Administration

https://www.cdc.gov/injection-safety/hcp/infection-control/?CDC_AAref_Val=https://www.cdc.gov/injectionsafety/fingerstick-devicesbgm.html

3. Self-Monitoring Blood Glucose Test Systems for Over-the-Counter Use

<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/self-monitoring-blood-glucose-test-systems-over-counter-use>

1. Use this device **ONLY** for the intended use described.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. This device does **NOT** serve as a cure for any symptoms or diseases.
The data measured is for reference only.

5. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
6. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.), may cause damaging static discharges that may cause erroneous results.
7. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
8. Do **NOT** under any circumstances use the device on neonates or infants.
9. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
10. Proper maintenance and periodically control solution test are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact your local customer service or place of purchase for help.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

BEFORE YOU BEGIN	07
Important Information	07
Intended Use	08
Test Principle	08
Contents of System	09
Meter Overview	10
Display Screen	11
Test Strip	12
Setting the Meter	13
THE FOUR MEASURING MODES	17
Control Solution Testing	18
Performing a Control Solution Test	18
TESTING WITH BLOOD SAMPLE	21
Preparing the Puncture Site	21
Alternative Site Testing	23
Performing a Blood Glucose Test	25
METER MEMORY	28
Reviewing Test Results	28
Reviewing Blood Glucose Day Average Results	29
Data Transmission Via Bluetooth	29
MAINTENANCE	31
Battery	31
Caring for Your Meter	32
Disinfecting Procedures	33
Caring for Your Test Strips	35
Important Control Solution Information	36
SYSTEM TROUBLESHOOTING	37
Result Readings	37
Error Messages	38
Troubleshooting	40
DETAILED INFORMATION	41
SPECIFICATIONS	42
FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT	43

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have any symptoms of illness, first repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycaemic-hyperosmolar state, with or without ketosis. Please consult your healthcare professional before use.
- Critically ill patients should not be tested with a blood glucose monitoring system.

Intended Use

The **FORA TN'G Voice** Blood Glucose Monitoring System is intended for use in the quantitative measurement of glucose in fresh capillary whole blood from the finger and the following alternative sites: palm, forearm or upper-arm. It is intended for use by people with diabetes mellitus at home as an aid in monitoring the effectiveness of diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and is not intended for use on neonates. It is intended to be used by a single person and should not be shared.

The alternative site testing in the **FORA TN'G Voice** Blood Glucose Monitoring System can be used only during steady-state blood glucose conditions.

This meter has some speaking functions but is not intended for use by the visually impaired.


Test Principle

Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

This system uses enzyme of Glucose Dehydrogenase (GDH) from *E. coli*.

Contents of System

The **FORA TN'G Voice** system kit includes:


- 1 Meter
- 2 Owner's Manual
- 3 Protective Wallet
- 4 Quick Start User Guide / Warranty Card
- 5 Daily Log Book
- 6 Two AAA alkaline batteries

Test strip, control solution, lancing device or sterile lancets may be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

NOTE:

If any items are missing from your kit or opened prior to use, or to obtain test strips and control solutions, please contact your place of purchase or call Customer Service for assistance.

Meter Overview

1 Test Slot with Strip Indicator

Insert test strip here to turn the meter on for testing.

2 Test Strip Ejector

Eject the used strip by pushing up this button.

3 Main Button

Enter the meter memory and mute a reminder alarm.

4 Display Screen

5 Bluetooth Switch

Activating the Bluetooth function.


6 SET Button

Enter and confirm the meter settings.

7 Bluetooth Indicator

Download test results with a Bluetooth connection.

Display Screen

- | | |
|-------------------------|-----------------------|
| 1 Volume Symbol | 9 Time |
| 2 Test Result | 10 Day Average |
| 3 Ketone Warning | 11 Alarm Symbol |
| 4 Warning Symbol | 12 Measurement Unit |
| 5 Testing Mode Symbol | 13 Low Battery Symbol |
| 6 Memory Symbol | 14 Test Strip Symbol |
| 7 Control Solution Mode | 15 Blood Drop Symbol |
| 8 Date | |

Test Strip

Absorbent Hole

Collect a drop of blood here. The blood will be automatically absorbed.

Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

ATTENTION:

The front side of test strip should face up when inserting test strip.

Test results might be wrong if the contact bar is not fully inserted into the test slot.

NOTE:

The **FORA TN'G Voice** Blood Glucose Monitoring System should only be used with **FORA TN'G Voice** Test Strips. Using other test strips with this monitoring system can produce inaccurate results.

Setting the Meter

Before using your meter for the first time or if you change the meter battery, you should check and update these settings. Make sure you complete the steps below and have your desired settings saved.

Entering the Setting Mode

Start with the meter off (no test strip inserted). **Press SET.**

1. Setting the date

With the year flashing, press Main button until the correct year appears. Press SET.

With the month flashing, press Main button until the correct month appears. Press SET.



With the day flashing, press Main button until the correct day appears. Press SET.

2. Setting the time format

Press Main button to select the desired time format — 12h or 24h. Press SET.

3. Setting the time

With the hour flashing, press Main button until the correct hour appears. Press SET.


With the minute flashing, press Main button until the correct minute appears. Press SET.

4. Deleting the memory

With “dEL” and a flashing “” on the display, press Main button and select “no” to keep the results in memory then press SET to skip.

To delete all the results, press Main button and “yes” and “” are displayed on the meter, press SET to delete the memory.


5. Setting the speaking volume

There are seven (7) speaking volume options to choose from. Press Main button until the desired speaking volume appears. To confirm your selection, press SET.


Volume 0 indicates that the speaking function is turned off, and “” will not display during testing. Volume 1 to 7 indicates speaking volume from low to

high, and “” will be displayed throughout the testing.


6. Choosing a language

Press the Main button to select L1/ L2. The default language for the meter is L1, which is English. To confirm your selection, press SET.

7. Setting the reminder alarm

You may set up any or all of the reminder alarms (1-4). The meter displays “On” or “OFF” and “”; press Main button to turn on or turn off to set the first reminder alarm.


Press Main button to select “On”, and then press SET to set the hour. When the hour is flashing, press Main button to increase an hour. Press SET to confirm and go to minutes, and then press Main button to increase one minute. Hold Main button longer to add faster. Press SET to confirm and go to the next alarm setting.

If you do not want to set an alarm, press SET to skip this step.

If you want to turn off an alarm, select the alarm number by pressing SET in the setting mode, and then press Main button to change displayed message from “ON” to “OFF”.

At the pre-set alarm time, the meter will beep and automatically turn on. You can press Main button to mute the alarm and insert a test strip to begin testing. If you do not press Main button, the meter will beep for 2 minutes and then switch off. If you do not want to test at this time, press Main button to switch off the meter.

8. Setting the Auto-Send

Press Main button to select the auto-send On or OFF. Press SET.

NOTE: This function is referring to the Bluetooth data transmission. If “On” is selected, your result will be transmitted automatically right after the test.

Congratulations! You have completed all settings!

NOTE:


- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

THE FOUR MEASURING MODES

The meter provides you with four modes for measuring, General, AC, PC and QC.

MODES	USE WHEN
General (displays as “Gen”)	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	performing a control solution test



You can switch between each mode by:



1. Start with the meter switched off.
Insert a test strip to turn on the meter.
The screen will display a flashing “” and “Gen”.

Voice Guide

Please apply blood onto the strip.

2. Press Main button to switch between General, AC, PC and QC mode.

Control Solution Testing


Our Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.

Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

Do a control solution test when:

- you first receive the meter,
- at least once a week to routinely check the meter and test strips,
- you open a new vial of test strips,
- you suspect the meter or test strips are not working properly,
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- practicing the testing process, or
- you have dropped or think you may have damaged the meter.

Performing a Control Solution Test



1. Insert the test strip to turn on the meter

Insert the test strip into the meter. Wait for the meter to display “” and “”.



Voice Guide

Thank you for using this product. Please relax during measurement. Please apply blood onto the test strip.



2. Press Main button to mark this test as a control solution test

With “**QC**” displayed, the meter will not store your test result in memory. If you press the Main button again, the “**QC**” will disappear and this test is no longer a control solution test.


WARNING:

When doing the control solution test, you have to mark it so that the test result will **NOT** mix with the blood glucose **TEST RESULTS** stored in the memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

3. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

4. Read and compare the result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

(112 mg/dL = 6.2 mmol/L; 100–165 mg/dL = 5.6–9.2 mmol/L)

Voice Guide

The blood glucose is (number) milligram per deciliter.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Go to www.foracare.com or call customer service for help.

NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **Meter Maintenance** section for important information about your control solutions.

TESTING WITH BLOOD SAMPLE

Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.
- Wash and dry your hands thoroughly before and after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by that person.

Sharing the lancing device and lancets may increase the risk of spreading infectious diseases. Lancing device must not be used on more than one person.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

Blood from an alternate site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.


Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either at fingertips or another body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting up the lancing device.

● Fingertip testing

Press the lancing device tip firmly against the lower side of your fingertip. Press the release button to prick your finger, you will hear a click which indicates that the puncture is complete.


● Blood from sites other than the fingertip

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm or hand, avoid lancing the areas with obvious veins because of excessive bleeding.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing

Important:

There are limitations with AST (Alternative Site Testing).

- **AST results should not be used for CGM calibration.**
- **AST results should not be used for insulin dosing calculations.**

Please consult your health care professional before you perform AST.

What is AST?

Alternative site testing (AST) means that people use parts of the body other than the fingertips to check their blood glucose levels. This system provides you to test on the palm, the forearm or the upper arm with the equivalent results to fingertip testing.

What is the advantage?

Your fingertips feel more pain because they are full of nerve endings (receptors). At other body sites, where nerve endings are not so condensed, you will not feel as much pain as you would pricking your fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, take a blood sample from your finger only.


We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not match the way you feel.
- Your routine glucose results often fluctuate.

Performing a Blood Glucose Test

1. Insert the test strip to turn on the meter

Wait for the meter to display “” and “”.


Voice Guide

Thank you for using this product. Please relax during measurement. Please apply blood onto the test strip.


2. Select the appropriate measuring mode by pressing Main button

For selecting the measurement mode, please refer to the “THE FOUR MEASURING MODES” section. (page 17)

3. Obtaining a blood sample

Use the pre-set lancing device to puncture the desired site. Wipe off the first appeared drop of blood with a clean cotton swab. The size of the drop should be at least as big as  (actual size), which is **0.5** microliter (μL) of volume. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

4. Bring strip to sample

While holding the meter in an angled position with the strip facing down, gently bring the strip to the sample so that the absorbent hole can sip in the blood. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not bring the strip to the blood sample within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

5. Read your result

The result of your blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.

(100 mg/dL = 5.6 mmol/L)

Voice Guide

The blood glucose is (number) milligrams per deciliter.

6. Eject the used test strip

Eject the test strip by pushing the eject button on the side. Use a sharps container to dispose the used test strips. The meter will switch itself off automatically.

Always follow the instructions in the lancing device insert when removing the lancet.

Take the lancet out carefully. Place the disk on a hard surface and push the exposed tip into the protective disk.

When removing the lancet, follow the instructions from the lancing device insert.


WARNING:

The used lancet and test strip may be biohazardous . Please discard them properly.

METER MEMORY

The meter stores the 450 most recent blood glucose test results along with respective dates and times in its memory. To enter the meter memory, **start with the meter switched off.**


Reviewing Test Results

1. Press and release Main button


“” will appear on the display. Press Main button again, and the first reading you see is the last blood glucose result along with date, time and the measurement mode.

(300 mg/dL = 16.7 mmol/L)

2. Press the Main button to recall the test results stored in the meter each time you press. **After the last test result, press Main button again and the meter will be turned off.**



(200 mg/dL = 11.1 mmol/L)

Voice Guide


The blood glucose is (number) milligrams per deciliter.

Reviewing Blood Glucose Day Average Results

1. Press and release Main button.

When “**M**” appears on the display, keep pressing Main button for 3 seconds until the flashing “**DAY AVG**” appears. Release Main button and then your 7-day average result measured in general mode will appear on the display.
(100 mg/dL = 5.6 mmol/L)

2. Press Main button to review 14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode in the order of Gen, AC, then PC.

3. Exit the meter memory. Keep pressing the Main button and the meter will turn off after displaying the last test result.

NOTE:

- Any time you wish to exit the memory, hold the Main button for 5 seconds or leave the meter idle for 3 minutes. The meter will turn off automatically.
- Control solution results are **NOT** included in the day average.
- If there are no records in memory, “—” displays when you recall the test results or review the average results.

Data Transmission Via Bluetooth

You can transmit your data from the meter to your devices (e.g. smart phone, tablet, PC...) via Bluetooth. Please contact your local customer service or place of purchase for assistance. Please note that you must

complete the pairing between the meter and Bluetooth receiver before transmitting data.

Pairing with your mobile device

1. Press and hold Bluetooth switch on the meter until Bluetooth indicator starts to flash.
2. On your mobile device, turn on the Bluetooth function and follow the instructions on your APP to pair the device (e.g., scan for meter to add into APP).
3. After successfully pairing the APP to you meter, the Bluetooth function on the meter will then start to transmit your data to your APP.

Bluetooth indicator on the meter

BLUETOOTH INDICATOR	STATUS
Flashing Blue	The Bluetooth function is enabled and waiting for connection.
Solid Blue	The Bluetooth connection is established.

WARNING:

- While the meter is in transmission mode, it will be unable to perform a test.
- Make sure your iOS (6 or above) or Android system (4.3 or above) device has its Bluetooth enabled and your meter is within receiving range before you transmit data from your meter.
- A compatibility issue between your mobile device and the meter might occur due to the way Bluetooth functionality is implemented by the various mobile device manufacturers.

MAINTENANCE

Battery

Your meter comes with two 1.5V AAA size alkaline batteries.

Low Battery Signal

The meter will display one of the messages below to alert you when the meter power is getting low.

The “” appears with **E-b, Error and low:**

The power is not enough to do a test. Please change the batteries immediately.

Replacing the Battery

To replace the batteries, make sure that the meter is turned off.

1. Press the edge of the battery cover and lift it up to remove.
2. Remove the old batteries and replace with two 1.5V AAA size alkaline batteries.
3. Close the battery cover. If the batteries are inserted correctly, you will hear a “beep” afterwards.

NOTE:

- Replacing the batteries does not affect the test results stored in the memory.

- As with all small batteries, these batteries should be kept away from children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose of the batteries according to your local environmental regulations.

Caring for Your Meter

Wash and dry your hands thoroughly before use or handling of the meter and strips.

What is Cleaning and Disinfection?

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter?

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter?

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt, blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. **Do NOT use organic solvents to clean the meter.**



We recommend using the product listed below to clean and disinfect your meter. The recommended product has been shown to be safe when used on the **FORA TN'G Voice** blood glucose monitoring system.

► **Micro-Kill+™ (Micro-Kill Plus™) by Medline**
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

- 1.** Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

- 2.** Wipe all exterior surface of the meter including the display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot, be careful not to allow excess liquid to get inside. Keep the meter moist with disinfection solution contained in the wipe for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe. Use two or more wipes if necessary.

- 3.** Remove the wipe. Allow the meter surface to dry completely.
- 4.** Discard the used wipes and never reuse them.

Wash your hands thoroughly with soap and water after handling the blood glucose monitoring system, lancing device and test strips to avoid contamination.

Each cleaning and disinfection cycle includes a pre-cleaning step with one wipe and a disinfection step with a second wipe. Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-888-307-8188 for assistance.

This device has been validated to withstand 260 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 1 cleaning and disinfection cycle per meter per week for 5 years, the expected life of device.

Stop using the meter if you see any signs of deterioration such as the LCD display is cracked or cloudy, the buttons no longer function, or the outer casing is cracked. Please contact customer service at 1-888-307-8188 for a replacement meter if any of the signs of deterioration are noticed.

The meter should be replaced after the validated number of disinfection cycles or the warranty period, whichever comes first.

NOTE :

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot, the battery compartment or the eject button.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

Meter Storage

- Storage conditions: -4 °F to 140 °F (-20 °C to 60 °C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and heavy impact.
- Avoid direct sunlight and high humidity.

Cleaning and Disinfection of Your Lancing Device

For home users, the lancing device is reusable. Please disinfect your lancing device regularly.

- To clean the lancing device, use a disinfecting wipe to clean the lancing device.
- Disinfect the cap by placing it in 70% alcohol for 10 minutes and then allow it to air dry.
- Do **NOT** place the device in a dishwasher or use detergents.

Caring for Your Test Strips

- Storage conditions: 35.6 °F to 89.6 °F (2 °C to 32 °C), below 85% relative humidity. Do **NOT** freeze.
- Store your test strips in their original vial only. Do not transfer to another container.
- Store test strip packages in a cool dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiration date. This may cause inaccurate results.

- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only our control solutions with your meter.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test be done at room temperature 68°F to 77°F (20°C to 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 35.6°F to 86°F (2°C to 30°C). Do **NOT** freeze.

SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair yourself and never try to disassemble the meter under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
	Appears when your blood glucose test result is below the lower measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
	Appears when your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
	Appears when your blood glucose test result is higher than the upper limit of measurement, which is 600 mg/dL (33.3 mmol/L).

Error Messages

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	WHAT TO DO
	Only Error Beep	Appears when the batteries cannot provide enough power for a test.	Replace the batteries immediately.
	Only Error Beep	Appears when a used test strip is inserted.	Repeat with a new test strip.
	Only Error Beep	Appears when ambient temperature is below system operation range.	System operation range is 50°F to 104°F (10°C to 40°C). Repeat the test after the meter and test strip are in the above temperature range.
		Appears when ambient temperature is above system operation range.	
 	Only Error Beep	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.

MESSAGE	WHAT DOES METER SAY?	WHAT IT MEANS	WHAT TO DO
	Only Error Beep	Problem with the meter.	Repeat the test with a new test strip. If the meter still does not work, please contact the customer service for assistance.
			
			
			

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the battery.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer services.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 68°F to 77°F (20°C to 25°C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

DETAILED INFORMATION

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.

Please consult your doctor to determine a target range that works best for you.

SPECIFICATIONS

Model No.: FORA TN'G Voice

Dimension & Weight: 104.4 x 58 x 23 mm, 71.2 g

Power Source: Two 1.5V AAA alkaline batteries

Display: LCD

Memory: 450 measurement results with respective date and time

External Output: Bluetooth

Expected Service Life: 5 years

Auto electrode insertion detection

Auto reaction time count-down

Auto switch-off after 3 minutes without action

Temperature warning

Operating Condition:

50°F to 104°F (10°C to 40°C), below 85% R.H. (non-condensing)

Meter Storage/Transportation Conditions:

-4°F to 140°F (-20°C to 60°C), below 95% R.H.

Strip Storage/Transportation Conditions:

35.6°F to 89.6°F (2°C to 32°C), below 85% R.H.

Measurement Units: mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

Federal Communications Commission (FCC) Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including

interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

- 1) This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
- 2) This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

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Distributed by Foracare, Inc. / Distribuido por Foracare, Inc.