Blood Glucose Monitoring System G20

OWNER'S MANUAL

No code setting required

Questions?

Toll Free: 1-866-469-2632

Hours of operation: 24 hours, 7 days/week

MADE WITH PRIDE AND CARE FOR H-E-B, SAN ANTONIO, TX 78204

MADE IN TAIWAN

GUARANTEE

We believe the high quality of this H-E-B product makes it an outstanding value. We hope you'll agree. If not, we'll cheerfully refund your money. Thanks for shopping with us.

H-E-B Recommended by
Pharmacy

**FORA G20
BLOOD GLUCOSE
MONITORING SYSTEM**

English

Version 6.0 2025/04

541-4251100-001

Thank you for choosing

FORA G20

BLOOD GLUCOSE MONITORING SYSTEM

Dear **FORA G20** System Owner:

This manual contains important information you must know about the system. Please read it thoroughly and carefully.

The greatest feature of this system is its strip-ejection function, which helps you remove blood filled strip without touching it with your hand. This may prevent any cross-contamination of potential biohazard materials.

Another unique feature is its alarm function, which alerts you to perform blood glucose test at the designated time. This thoughtful design helps you routinely check your daily blood glucose on time.

Another unique feature of this system is its no code function, which means that there is no need to calibrate your meter. You just simply insert the test strip into the meter and begin the test.

For other questions regarding this system, please contact our customer service team at 1-888-307-8188, 1-866-469-2632. You can also visit www.foracare.com for helpful information.

IMPORTANT SAFETY INSTRUCTIONS

READ THIS BEFORE USING

The following basic safety precautions should always be taken.

- 1.Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
- 2.Use the device only for the intended use described in this manual.
- 3.Do not use accessories which are not supplied by the manufacturer.
- 4.Do not let the device or its cord come into contact with surfaces that are hot to the touch.
- 5.Do not use the device near aerosol sprays, or where oxygen is being administered.
- 6.Do not use device if it is not working properly, or if it has suffered any damages.
- 7.Before using product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

IMPORTANT INFORMATION	06
ALTERNATIVE SITE TESTING (AST)	07
INTRODUCTION OF THE SYSTEM	09
Intended Use	09
Principle of Measurement	09
Contents of the System	10
Appearance and Key Function of the Meter	11
LCD Screen	12
Test Strip	13
PREPARATIONS BEFORE USE	14
Battery Replacement	14
Setting the Meter and Deleting the Memory	16
Alarm Function	24
BEFORE TESTING	27
Checking with FORA Control Solutions	27
Important Control Solution Information	28
Performing a Control Solution Test	29
TESTING YOUR BLOOD	31
Test Procedure	32
Expected Test Results	38
COMPARING METER AND LABORATORY RESULTS	39
USING THE METER MEMORY	40
VIEWING RESULTS ON A PERSONAL COMPUTER	43

TAKE CARE OF YOUR METER AND STRIP	44
Cleaning	44
Storage	44
TROUBLE SHOOTING GUIDE	46
Special Message	47
Error Message	48
Problem in Operation	49
SPECIFICATIONS	51

IMPORTANT INFORMATION

- Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If you get your blood glucose results lower or higher than usual, and do not have symptoms, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Apply only capillary whole blood sample to test your blood glucose. Applying other substances will cause wrong results.
- If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate low results may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.
- Please refer to your test strip package insert for additional important information.


ALTERNATIVE SITE TESTING (AST)

Important:

There are limitations for doing AST. Please consult your healthcare professional before you do AST.

What is AST?

Alternative site testing (AST) means that people use parts of the body other than fingertips to check their blood glucose levels. This system provides you to test on the palm, the forearm, the upper arm, the calf, and the thigh with the equivalent results to fingertip testing.

What's the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, since nerve endings were not so condensed, you will not feel as much pain as at the fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at fingertip reflects these changes faster than capillary blood at other sites. Therefore when testing blood glucose during or immediately after meal, physical exercise, or any other event, **take blood sample from your finger only.**

We strongly recommend you do AST **ONLY** in the following intervals:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unawareness of hypoglycemia.
- Your AST results do not match the way you feel.
- You are testing for hyperglycemia.
- Your routine glucose results are often fluctuating.
- You are pregnant.

How to increase the accuracy?

Stimulating blood perfusion by rubbing the puncture site prior to blood extraction has a significant influence on the glucose value obtained. Blood from the site without rubbing exhibits a measurably different glucose concentration than blood from the finger.

Please follow suggestions below before getting a drop of blood:

- Rub the puncture site about 20 seconds before penetration.
- Use a clear cap (included in the kit) instead while setting the lancing device.

INTRODUCTION OF THE SYSTEM

Intended Use

The system is intended for use outside the body (in vitro diagnostic use). It should be used only for testing glucose (sugar) and only with fresh capillary whole blood samples taken from the finger and the alternative sites including the palm, the forearm, the upper arm, the calf and the thigh. The system is intended for use in the home and in clinical settings. It should not be used for the diagnosis of diabetes or for the testing of newborns.

AST in this system can be used only during steady-state blood glucose conditions described in the section of “About AST”.

Principle of Measurement

The test is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

Contents of the System

The two main parts are included in this system: a blood glucose meter and test strips. These products have been designed, tested, and proven to work together as a complete system to measure the concentration of your blood glucose precisely.

Always use the FORA brand of test strips and a control solution with your FORA glucose meter.

Your system includes:


- | | | |
|--------------------------------|-------------------------|-----------------|
| 1 1 meter | 2 Owner's manual | 3 Warranty card |
| 4 Quick start user guide | 5 Daily log book | |
| 6 Protective wallet | 7 1 x CR2032 3V battery | |
| 8 10 blood glucose test strips | | |

NOTE

Please make sure that all products listed above are contained and sealed in the package before using this system. If you find any imperfection in our products, please return the whole system to the place of purchase.

Appearance and Key Function of the Meter

1. LCD DISPLAY

Guides you through the test using symbols and simple messages.

2. MAIN BUTTON

Located in front of the meter with “M” on it, is used to turn on the meter, enter the memory or control steps of setting.

3. TEST SLOT

Is where you insert the test strip. The meter will turn on automatically after insertion.

4. STRIP-EJECTION BUTTON

Is where the used strip will be ejected after you push up the button.

5. DATA PORT

Located at the side, is for cable connection.


6. C BUTTON

Located at the side, is for control solution test.

7. SET BUTTON

Located in the battery compartment, is used to set up the meter.

LCD Screen

BLOOD DROP SYMBOL

Flashes when it is ready to apply the sample.

TEST STRIP SYMBOL

Appears when the meter is turned on.

LOW BATTERY SYMBOL

Appears when the battery power is low

TEST RESULT AREA

MEASUREMENT UNIT

MEMORY SYMBOL

Appears when you review the memory.

DAY AVERAGE

Indicates that the displayed test result is an average.

ALARM FUNCTION INDICATOR

Appears when alarms are set.
Four alarms can be set in this meter.

KETONE WARNING

Appears when the test result is equal or higher than 240mg/dL(13.3mmol/L)

CTL SYMBOL

Appears when doing a control test and indicates that the result won't be stored in the memory.

DATE



TIME

Test Strip

The front side of the test strip should face up when inserting the test strip.

ATTENTION:

The device will not turn on if strip is inserted halfway or incorrectly.

NOTE

The **FORA G20** meter should only be used with **FORA G20** Test Strips. Using other test strips with this meter can produce inaccurate results.


PREPARATIONS BEFORE USE

Battery Replacement

Your meter comes with one 3V CR2032 Lithium battery. The meter will alert you when the power is low by displaying two different messages:

1. When  symbol is displayed on the main screen: the meter is functional and the result remains accurate, but battery should be replaced.


2. With  symbol, low and E-b symbols on the displays: the battery can not provide enough power to do a test. You must replace the battery immediately.

NOTE

- Replacing the battery within 5 minutes does not affect the settings. Therefore, the settings do not need to be updated when the battery is changed.
- Batteries might leak chemicals if not used for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more)
- Do not use the used battery. Use only 3V CR2032 new battery.


To replace the battery, make sure that the meter is turned off.

STEP1

STEP1


Press the buckle on battery cover and lift up to remove cover.

STEP2

STEP2

Remove the old battery and replace with one 3V CR2032 Lithium battery.

STEP3

STEP3

Close the battery cover.


WARNING

As with all small batteries, the batteries should be kept away from small children who still put things in their mouths. If they are swallowed, promptly see a doctor for help.

Setting the Meter and Deleting the Memory


Your meter comes with the time, date, unit of measurement, memory deletion and alarm function setting. If you need to reset these parameters, please follow the steps below.

Start with the meter off. Then press the set button located in the battery compartment. The meter is now in the setting mode.

STEP 1. Set the Year

After pressing the set button, the year, a blinking number, will be shown on the screen. Press and release the M button to plus one year. You can also keep pushing down the M button to proceed faster. When the correct year is displayed, press the set button and then a flashing number, which stands for the month, is shown.

STEP 2. Set the Month

Press and release the M button until you see the correct month. To move faster, keep pushing the M button down. When the desired month is displayed, press the set button and then the day will flash.

STEP 3. Set the Day

Press and release the M button until you see the correct day. To move faster, keep pushing the M button down. When the desired month is displayed, press the set button and then 12h or 24h will be shown on the screen.


STEP 4. Select the 12h or 24h

Press and release the M button to select 12h or 24h. It differs in time presentation format. If you select 12h, AM and PM will be shown together with time. An example is shown below.

If you select 24h, format of time will be from 00:00 to 23:59 without AM and PM. An example is shown below.

After finishing setting 12h or 24h, press the set button and then the hour will flash on the screen.

STEP 5. Set the Hour

Press and release the M button until you see the correct hour on the screen. To move faster, keep pushing the M button down. When the desired hour is displayed on the screen, press the set button and then the minute will flash.


STEP 6. Set the Minutes

Press and release the M button to advance one minute. To move faster, keep pushing the M button down. When the desired minute is displayed, press the set button to go the next step for "Delete Memory".

STEP 7. Delete Memory

"dEL", "[M]" and flashing Yes/no is displayed on the screen. If you do **NOT** want to delete memory, press the M button to select "no" and then press set button to go to step 9. If you'd like to delete ALL memory, press M button to select "yes". Then **press set button to delete All memory**. "OK" is displayed in the meter, which means that all data stored is deleted. Go to step 9.

NOTE


- The time, date, memory deletion function and alarm function can **ONLY** be changed in the setting mode. Therefore, during glucose testing, those parameters cannot be changed.
- The meter cannot automatically update daylight saving time. You have to manually adjust the time in the meter according to the procedures. Please make sure the meter matches the time zone and adjust the daylight saving time data if required.
- Your meter displays 7-, 14-, 21-, 28-, 60- and 90-day averages which you can access from the meter memory. These averages are calculated from results obtained during the 7-, 14-, 21-, 28-, 60- and 90-day preceding the current date and time settings. When the date and time are changed, the 7-, 14-, 21-, 28-, 60- and 90-day averages may change.
- While the meter is in the setting mode, it will turn off automatically without any action in three minute.

STEP 8. Select and Set Alarm Function


The meter provides four alarms. You are able to set four different alarming times a day. "On" or "OFF" and "1" are displayed on the screen firstly. Press the M button to turn on or turn off Alarm 1.

► Alarm 1 setting

- If you press the M button to select "OFF", then press the set button to go to set Alarm 2.


- If you press the M button to select "On", then press set button to set hour. While the hour is flashing, press the M button to plus an hour, and then press the set button. Blinking minute is displayed on the screen. Press the M button to advance one minute. To move faster, keep pushing the M button down. When the desired minute is displayed, press the set button to proceed to the next alarm setting.


► Alarm 2 setting

- If you press the M button to select “OFF”, then press the set button to go to set Alarm 3.


- If you press the M button to select “On”, then press set button to set hour. Set hour and minute according to the instruction described in Alarm 1. After finished, press the set button to go to set Alarm 3.


► Alarm 3 setting

- If you press the M button to select “OFF”, then press the set button to go to set Alarm 4.


- If you press the M button to select “On”, then press set button to set hour. Set hour and minute according to the instruction described in Alarm 1. After finished, press the set button to go to set Alarm 4.


► Alarm 4 setting

- If you press the M button to select “OFF”, press the set button to turn off the meter.

- If you press the M button to select “On”, then press set button to set hour and minute according to the instruction described in Alarm 1. After finished setting Alarm 4, press set button to turn off the meter.

Alarm Function


The meter provides for alarms. The alarm function alerts you to perform your blood glucose regularly. If you need to set alarm function, please refer to page20 " Select and Set Alarm Function" for information.

If you have turned on the alarm function, the meter will auto-start in the time point you set with short beeps for two minutes.


An example if set Alarm 1 at 10:00AM.

The meter automatically turns on at 10:50 AM everyday with beeping for two minutes, blinking "  " and . You can directly insert a new strip to perform blood glucose test now.


An example if Alarm 1 at 10:00AM and Alarm 2 at 2:00PM.

The meter automatically turns on at 10:00 AM first and next turns on at 2:00 PM. You can directly insert a new strip to perform a test.

If you don't want to perform a test after the meter turns on:

- Press the M button. "OFF" displays and then the meter will turn off.

- Or you can leave the meter without pressing any button for three minutes. The meter will turn off automatically.

NOTE

- If the alarm function is set, the alarm symbol will appear on the screen when the meter is off.
- Each alarm setting function is independent, you can set alarm function separately and the setting will not affect other alarm function.

BEFORE TESTING

Checking with FORA Control Solutions

FORA control solutions contain a known amount of glucose that reacts with test strips. By comparing your control solution test results with the expected range printed on the test strip vial label, it is able to check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you do this simple check routinely to make sure you get accurate results.

How often should the control solution test be performed?

- When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- For routinely check the meter and test strips, perform a single test for each level of control solution at least once a week.

When should the control solution test be performed?

- When you first get your Glucose Meter.
- When you begin using a new vial of test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are not accurate.
- When your test strips are exposed to extreme environmental conditions (See **Storage** section of this manual).
- When you want to practice running the test.
- If you drop the meter.

Important Control Solution Information

- Use only FORA control solutions.
- Check the expiration date on the control solution vial. Do not use if expired.
- Control solution, meter, and test strips should come to room temperature 68 to 77°F (20 to 25°C) before testing.
- Shake the vial, discard the first drop of control solution, and wipe off the dispenser tip to ensure a good sample and an accurate result.
- Use only for 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.
- Store the control solution tightly closed at temperatures 36 to 86°F (2 to 30°C). Do not freeze.

NOTE


The control solution range is printed on FORA test strip vial. Only use FORA control solution. It is used to test meter and test strip performance. It is not recommended range for your blood glucose level.

Performing a Control Solution Test

TAKE A TEST STRIP OUT WITH CLEAN AND DRY HANDS FIRST.


STEP 1. Insert the Test Strip

Insert a test strip with contact bars end first and facing up, into the test slot. The meter turns on automatically and displays the followings in sequence: “CHK” and “” → , flashing “” with date and time.

STEP 2. Press the C button

While the “” symbol is flashing, press the C button and you will see the “CTL” appear, which means that the meter is in the “Control Solution Testing Mode”. Under this mode, the test result of the control solution is not stored in the memory. If you decide not to perform a control solution test, press the C button again and the “CTL” sign will disappear.


WARNING

- Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.
- Every time you perform a control solution test, you must enter into the “CTL” test mode so that the test result will not be stored in the meter memory. Failure to do so will confuse the blood glucose test result with the control solution test result in memory.

STEP 3. Obtain Control Solution


Shake the control solution vial well. Remove the cap from the control solution bottle. Place cap on flat surface. Squeeze the vial, discard the first drop, and wipe off the dispenser tip to prevent contamination. Squeeze the vial again to get another drop and **apply the drop to the top of cap**.


STEP 4. Apply Control Solution

While holding the monitor, move the absorbent hole of the test strip to touch the drop of control solution. Then the drop will be automatically drawn into the test strip. Make sure the confirmation window is completely filled. The meter begins to count down.

To avoid contaminating the control solution with the content of the test strip, you have to place a drop of control solution on a clean surface. Do not directly apply control solution into a strip.

STEP 5. Read and Compare the Result


Out-of-range results

If test results falls outside the range printed on the test strip vial, check the section "Problem in Operation" in troubleshooting guide and repeat the test. If you continue to get out-of-range results, it means that the system may not be working properly. **Do NOT** test your blood. Contact the Customer Care Line at 1-888-307-8188, 1-866-469-2632

TESTING YOUR BLOOD

Be sure to read this section and the test strip package insert found in the test strip box carefully before testing. Make sure you have all items needed to test:

- A. Blood Glucose Meter
- B. Test Strip
- C. Lancing Device
- D. Sterile Lancet
- E. Clear Cap (For AST use)

Lancing device, lancets, and clear cap may not be included.

WARNING

To reduce the chance of infection:


- Never share lancets or the lancing device with others.
- Always use a new sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.

Test Procedure

Wash and dry your hands before starting.

Step 1 Setting up the Lancing Device

Pull off the cap of the lancing device.


Insert a lancet into the lancet holder and push down firmly until it is fully secured.


Twist the protective disk off the lancet.

Replace the cap by aligning the arrow on the cap with the release button.


Select the depth of penetration by turning the adjustable tip in either direction so that the arrow on the cap points to the desired depth.

Pull the cocking control back until it clicks.
You will see a color change inside the
release button when it is ready.

If it does not click, the device may have
been cocked when the lancet was inserted.


The lancing device is now ready for use. Set aside for later use.

STEP 2. Insert the Test Strip

Insert a test strip by putting the contact bars end first and facing up into the test slot. The meter turns on automatically and shows the following messages:

“CHK” and “ → ”, flashing “” with date and time.

STEP 3. Get a Drop of Blood

Select the puncture site either in finger or in other parts (AST). Clean the puncture site with 70% alcohol cotton and **let it air dry**.

- **Fingertip Testing**


Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.

• Blood from Sites Other Than the Fingertip

Please refer to the section of “About AST” for available punctured sites.

After penetration, gently massage the punctured area to obtain blood. The volume of blood sample must be at least 0.5 microliter (•actual size). But be attention **NOT to smear the blood sample**. It is recommended to discard the first drop of blood with a clean tissue paper or cotton and use the second drop of blood for testing.


NOTE

- Choose a different spot each time you test. Repeated punctures in the same spot may cause soreness and calluses.
- Before you decide to do AST, please consult your health professional firstly.
- Since the first drop of blood usually contains tissue fluid and serum, which may affect the test result, it is recommended to be discarded.


STEP 4. Apply Blood into the Test Strip

When “💧” is flashing on the screen, apply your blood to the absorbent hole of the test strip until the confirmation window is fully filled with blood. The meter then begins to count down automatically.

STEP 5. Obtain an Accurate Result in 7 Seconds

The result of your blood glucose test is shown after the meter counts to 0. This reading is automatically stored in the meter.


NOTE

- Do not push your finger (with blood on it) against the test strip or try to apply a smeared sample on the test strip.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should completely fill the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting, **NEVER** try to add more blood to the test strip. **Discard the test strip and retest with a new one.**
- If you have trouble filling the test strip, please contact the customer service for help.

STEP 6. Eject the Used Test Strip


After finishing the measurement, you can either take out the used strip by simply sliding up the Strip-Ejection button or remove the test strip directly with your hand. "OFF" is shown after the used strip is ejected and the meter will shut off automatically afterward.

STEP 7. Remove the Lancet

Always use caution when removing the lancet.

Take the lancet out carefully. Place the disk on a hard surface and push the exposed tip into the protective disk.

WARNING

The used lancet and the used test strip may be potentially biohazard. Please discard it carefully according to your local regulations.

Expected Test Results

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that keeping blood glucose levels close to normal can reduce the risk of diabetes complications by up to 60%*¹. The results you get with the **FORA G20** system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)

* Source: American Diabetes Association (2008). Clinical Practice Recommendations. Diabetes Care, 31 (Supplement 1): S1-108.

*¹: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Please work with your doctor to determine a target range that works best for you.

COMPARING METER AND LABORATORY RESULTS

The meter provides you with whole blood equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way (See test strip package insert for typical accuracy and precision data, and for important information on Limitations). To make an accurate comparison between meter and laboratory results, follow the guide-lines below.

Before you go to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- It is best to fast for at least eight hours before doing comparison tests.
- Take your meter with you to the lab.

While at the lab:

Make sure that the samples for both tests (the meter test and the lab test are taken and tested within 15 minutes of each other).

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication, or experienced stress^{*2}. In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test^{*3}. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (severe dehydration) may also cause a meter result to be different from a laboratory result.

References

*2: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

*3: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

USING THE METER MEMORY


Your meter stores the latest 450 results of your blood glucose with date and time in its memory. It also automatically calculates the averages of your blood glucose in the intervals of 7, 14, 21, 28, 60 and 90 days. You can review the results easily by the following steps.

Recall the Stored Test Results

STEP 1.

When the meter is off, press and release the M button.

The screen shows "M". Press the M button again, "01" appears first and then the latest glucose result along with date and time will be shown on the screen.


STEP 2.

Press the M button once by once to recall the test results stored in the meter consecutively.

STEP 3.

After the last test result, press the M button again and the meter will be turned off.


Read the Average of Blood Glucose Results


STEP 1.

When the meter is off, press and release the M button.

The screen shows “M”. Keep pressing the M button for 2-3 seconds, until blinking “DAY AVG” appears. Release the M button and then the 7-day average result will appear on the screen.

The 7-day average is calculated from the blood glucose results obtained during the last 7 days.

1. The **average** of glucose tests.
2. The average was calculated from the test results of the **last 7 days**.
3. **6 glucose tests** have been performed in the last 7 days.
4. You can interpret the figure as: The average of 6 glucose test in the last 7 days is 75 mg/dL

STEP 2.




Press the M button once by once

to review the 14-, 21-, 28-, 60- and 90- day average in order. Like the 7-day average, the 14- day average and the times of performing test in the past 14 days will be shown on the screen.

STEP 3.

After the 90-day average, press the M button again to resume and it will go to the most recent test result. Follow the steps in “Recall the Stored Test Results” to review individual result.

NOTE

- When pressing the M button to recall the stored test results, blinking "M" will be displayed on the main screen first. If you do not press the M button within 5 seconds, the blinking " " will appear. In the meantime, you can choose to insert a test strip to **start testing your blood** or to press the M button again to review the **stored test results**.
 
- When using the meter for the first time, "---" is displayed when you recall the test results or review the average result. It means that there is no test result in memory.
 
- Anytime when you want to exit the memory, keep pressing the M button for 3 second until "OFF" displays. The meter then automatically shut down.
 
- If no button is pressed within 3 minutes, the meter will show "OFF" and turn off automatically.
- The control solution results are **NOT** stored in the memory (please also go to page **26 WARNING** for information). The list of past results and the average result are for blood glucose results only.

VIEWING RESULTS ON A PERSONAL COMPUTER


Results in memory can be transmitted to the personal computer. Health Care Software System and an Interface Cable are required before installation. The software can be downloaded from the home page of FORA, Inc. The interface cable is an optional accessory. To learn more about Health Care Software System or to obtain an Interface Cable separately, please contact your Customer Care Line at 1-888-307-8188, 1-866-469-2632.

Step 1. Installing the Software

Install Health Care Software System on your Computer by following the instructions provided on Fora Care's website: www.foracare.com

Step 2. Connecting to a Personal Computer

Connect the interface cable to the serial port in your computer. With the meter turned off, connect the Interface Cable to the Data Port of the meter. "PC" will appear on the display, indicating that the meter is ready to transmit data.

Step 3. Data transmission

Follow the instructions provided in the software to transmit data. Results transmitted will include date and time. Remove the cable and the meter will automatically turn off.

NOTE

While the meter is connected to the PC, it is unable to perform a blood glucose test.

TAKE CARE OF YOUR METER AND STRIP

To avoid the meter and test strips getting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

Cleaning

1. To clean the meter exterior, wipe with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft and dry cloth. Do not flush with water.
2. Do not use organic solvents to clean the meter.

Storage


1. Meter Storage


- Storage condition: -4°F to 140°F (-20°C to 60°C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and strong impact.
- Avoid direct sunlight and humidity.

2. Strip Storage

- Storage condition: 39.2 to 104°F (4 to 40°C), below 85% relative humidity. Do not freeze.
- Store your test strips in their original vial only. Do not transfer to other container.


- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately replace the vial cap and close it tightly.
- Handle the test strip with clean and dry hands.

- Use each test strip immediately after removing it from the vial.
- Write the discard date (the date opened plus 90 days) on the vial label when you first open it. Discard remaining test strips on this date.

- Do not bend, cut, or alter test strips in any way.
- Keep the strip vial away from children as the cap and the test strip may be choking hazards. If swallowed, promptly see a doctor for help.

3. Control Solution Storage

- Storage condition: Store the control solution tightly closed at temperatures 36 to 86°F (2 to 30°C). Do not freeze.
- Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.

TROUBLE SHOOTING GUIDE

The “Error Messages” you will probably encounter when using this meter are listed as followings. If you get any problem in performing test with this meter, please refer to the following problem-solving guides. These guides help you to identify and solve certain problems, but not in all circumstances. Improper use may cause inaccurate results without showing any error message or symbol. In the event of a problem, refer to the information under action.

Never try to disassemble the meter in any circumstances. If you encounter any error messages not listed below or if you have followed the actions recommended below but the problem keeps unsolved, please Contact the Customer Care Line at 1-888-307-8188, 1-866-469-2632.

Special Message

Special symbols and messages appear together with your test result.

MESSAGE	WHAT IT MEANS
Lo	appears when your result is below 20 mg/dL (1.1 mmol/L).
KETONE?	is shown when your result is equal or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for Type 1 diabetes. Please seek medical assistance immediately.
H i	is displayed when your result is higher than 600 mg/dL (33.3 mmol/L).

Error Message

MESSAGE	WHAT IT MEANS	Action
 Error low	Appear when the battery can not provide enough power for a test.	Replace the battery immediately.
 Error	Appear when inserting a used test strip.	Test with a new test strip.
 Error low	Appear when environmental temperature is below system operation range: 50°F(10°C)	System operation range is 50°F to 104°F(10°C to 40°C). Repeat the test after the meter and test strip have reached the above temperature.
 Error high	Appear when environmental temperature is above system operation range: 104°F(40°C).	
 Error	Remove the strip after applying blood to the absorbent hole.	Re-test with a new test strip.
 Error Error	Problem with the meter.	Review the instructions and re-test with a new test strip. If the above steps do not work, please contact the dealer.

Problem in Operation

1. If the meter does not display a message after inserting a test strip:

PROBABLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the battery.
Battery incorrectly installed or absent.	Check that the battery is correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the contact bars end first and "▲" facing up.
Defective meter.	Please contact the customer service for assistance.

2. If the test does not start after applying the sample:

PROBABLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.
Defective test strip.	Repeat the test with a new test strip.
Did not apply sample until "▲" is flashing.	Repeat the test with a new test strip. Apply sample only when "▲" appears on the display.
Defective meter.	Please contact the customer service for assistance.

3. If the control solution test result is out of range.

PROBABLE CAUSE	WHAT TO DO
Error in performing the test.	Read the instruction thoroughly and repeat the test again.
Do not shake the control solution vial very well.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date or the discarded date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should come to room temperature (68°F to 77°F(20°C to 25°C) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Meter malfunction.	Please contact the customer service for assistance.

SPECIFICATIONS

Model No.: TD-4251

Dimension & Weight: 86.9mm(L) x 51mm(W) x 15.8mm(H), 42.2g

Power source: One CR2032

Memory: 450 measurement results with date and time

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Auto turn-off after 3 minutes without action

Operating condition: 50°F ~104°F(10°C~40°C), below 85% R.H.(noncondensing)

Meter storage/Transportation condition: -4°F~140°F(-20°C~60°C), below 95% R.H.

Measurement Units: Fix mg/dL

Measurement Range: 20~600mg/dL (1.1~33.3mmol/L)

The device has been certified to meet the electrical and safety requirements of: IEC 60601-1, EN 60601-1, IEC 61010-1, EN 61010-1, EN 61010-2-101, EN 60601-1-2, EN 61326.

Performance characteristics

Accuracy

Within ± 15 mg/dL (0.83 mmol/L) at glucose concentration < 75 mg/dL (4.2 mmol/L) and within $\pm 20\%$ at glucose concentration ≥ 75 mg/dL (4.2 mmol/L).

Precision

CVs (%) of intermediate precision and repeatability are less than 5%.

The product has been tested to meet the requirements of ISO 15197.

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