
FORA[®] D40g Plus

Blood Glucose Plus Blood Pressure Monitoring System
*Sistema de Monitoreo de Glucosa y Medidor de
Presión Arterial*

Owner's Manual
Manual del Propietario

Dear FORA D40g Plus System Owner:

Thank you for purchasing the **FORA D40g Plus** Blood Glucose plus Blood Pressure Monitoring System. This manual provides important information which helps you to operate this system smoothly. Before using this product, please read the following contents thoroughly and carefully.

This system is intended to be used by a single patient. According to the clinical studies from the American Diabetes Association, elevated blood pressure often accompanies adult diabetes patients. These studies also suggest that diabetes patients could reduce the risk of cardiovascular diseases by managing their blood glucose levels and blood pressure. Therefore, with easy operation of this **FORA D40g Plus** Blood Glucose plus Blood Pressure Monitoring System, you can easily monitor your blood glucose levels and blood pressure by yourself at any place, any time. In addition, this system can help you and your healthcare professionals to monitor and adjust your treatment plans, and keep your diabetes and blood pressure under control.

If you have other questions regarding this product, please contact the place of purchase.

IMPORTANT SAFETY PRECAUTIONS

READ BEFORE USE

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. “Letter to Manufacturers of Blood Glucose Monitoring Systems Listed With the FDA”

<https://public4.pagefreezer.com/browse/FDA/12-01-2022T02:57/>

<https://www.fda.gov/medical-devices/in-vitro-diagnostics/letter-manufacturers-blood-glucose-monitoring-systems-listed-fda>

2. “CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens”

<https://www.cdc.gov/injectionsafety/fingerstick-devicesbgm.html>

3. “Infection Prevention during Blood Glucose Monitoring and Insulin Administration”

<https://www.cdc.gov/injectionsafety/blood-glucose-monitoring.html>

4. “Self-Monitoring Blood Glucose Test Systems for Over-the-Counter Use”

<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/self-monitoring-blood-glucose-test-systems-over-counter-use>

1. Use this device **ONLY** for the intended use described.
2. Do **NOT** use accessories which are not specified by the manufacturer.
3. Do **NOT** use the device if it is not working properly or if it is damaged.
4. Do **NOT** use the equipment where aerosol sprays are being used or where oxygen is being administered.

5. Do **NOT** use under any circumstances on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your physician to have the results interpreted.
7. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed.
8. Keep the device and testing supplies away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. Keep the device and its flexible cord away from hot surfaces.
10. Do **NOT** apply the cuff to areas other than the place directed.
11. Use of this device in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
12. Do **NOT** use this device in close proximity to sources of strong electromagnetic radiation, as these may interfere with the correct operation.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

TABLE OF CONTENTS

BEFORE YOU BEGIN	08
Important Information	08
Intended Use	09
Test Principle	10
Contents of System	11
Meter Overview	12
Display Screen	14
Test Strip	15
SETTING THE METER	16
THE FOUR MEASURING MODES	18
Control Solution Testing	19
Performing a Control Solution Test	19
TESTING WITH BLOOD SAMPLE	21
Preparing the Lancing Device for Blood Testing	22
Preparing the Puncture Site	22
Performing a Blood Glucose Test	23
Alternative Site Testing	26
TESTING YOUR BLOOD PRESSURE	27
Before Measurement	27
Fitting the Cuff Correctly	27
Proper Measurement Position	29
Taking Measurements	29
METER MEMORY	32
Reviewing Test Results	32
Reviewing Day Average Results	33
DOWNLOADING RESULTS	34

MAINTENANCE	37
Battery	37
Caring for Your Meter	40
SYSTEM TROUBLESHOOTING	45
Result Readings	45
Error Messages	46
Troubleshooting	48
DETAILED INFORMATION	50
SPECIFICATIONS	52
FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT	54

BEFORE YOU BEGIN

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to incorrect results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients in shock. Please consult the healthcare professional before use.

Intended Use

FORA D40g Plus Blood Glucose plus Blood Pressure Monitoring System is intended to be used for the quantitative measurement of glucose (sugar) in fresh capillary whole blood samples drawn from the fingertips, palm, forearm, and upper arm. This system is intended to be used by a single person and should not be shared.

The system is intended for self testing outside the body (*in vitro* diagnostic use) at home. This system should not be used for the diagnosis of or screening for diabetes, nor for neonatal use.

Alternative site testing should be done only during steady – state times (when glucose is not changing rapidly).

This system is also intended to be used as a non-invasive measuring system of systolic and diastolic blood pressure and pulse rate of an adult individual at home. The blood pressure is measured by using an inflatable cuff wrapped around the arm.

The meter contains some speaking functions but it has not been validated for use by visually impaired users.

Test Principle


Your system measures the amount of sugar (glucose) in whole blood. The test is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

The blood pressure is measured non-invasively at the arm based on the Oscillometric method.

This device is **NOT** able to take measurements in the presence of common arrhythmia, such as atrial or ventricular premature beats or atrial fibrillation. It may produce reading errors.

Contents of System

Your new **FORA D40g Plus** system kit includes:

- | | |
|---|-----------------------------------|
| 1 Meter | 5 Daily Log Book |
| 2 Owner's Manual | 6 One Rechargeable Li-Ion Battery |
| 3 Protective Wallet | 7 Pressure Cuff |
| 4 Quick Start User Guide /
Warranty Card | 8 AC Power Adaptor |


Test strip, control solution, sterile lancet and lancing device with one clear cap may be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand. Please contact the place of purchase or call customer service at 1-888-307-8188 for availability.

NOTE:


If any items are missing from your kit or opened prior to use, please contact local customer services or place of purchase for assistance.

Meter Overview


Front


Back

Side


Pressure Cuff

- 1 Display Screen**
- 2 RF Button**
Assists you with wireless data transmission.
- 3 S Button**
To enter and confirm the meter settings.
- 4 M Button**
To enter the meter memory.
- 5 AVERAGE Button**
Assist you with blood pressure average mode.
- 6 Speaker**
- 7 On/Off Button**
To start a single non-invasive blood pressure measurement.
- 8 Power Indicator**
Indicates the power status of the device.
- 9 USB Indicator**
Indicates when an USB plug is inserted.
- 10 RF Indicator**
Indicates a wireless data transmission.
- 11 Test Slot**
Insert test strip here to turn the meter on for testing.
- 12 Battery Compartment**
- 13 Strip Ejector**
Eject the used strip by pushing up this button.
- 14 AC/DC Adapter Port**
Connect to a power supply.
*Use the AC adapter provided by ForaCare only.
- 15 USB Link Port**
- 16 Pressure Cuff**
- 17 Air Tube**
- 18 Air Plug**
Connect to air jack.
- 19 Air Jack**

Display Screen

- | | | | |
|---|---|----|--------------------------------|
| 1 | Date | 10 | Pulse Rate |
| 2 | Systolic Pressure Symbol | 11 | Day Average Result |
| 3 | Measuring Mode | 12 | Battery Symbol |
| 4 | Systolic Pressure Value | 13 | Memory Mode Symbol |
| 5 | Diastolic Pressure Symbol | 14 | Unit for Blood Pressure |
| 6 | Diastolic Pressure Value | 15 | Unit for Blood Glucose |
| 7 | Volume Symbol | 16 | Time |
| 8 | Ketone Warning | 17 | Test Strip Symbol |
| 9 | Heart Symbol- for IHB
Detection (Irregular Heart
Beat) | 18 | Blood Drop Symbol |

Test Strip

Absorbent Hole

Apply a drop of blood here. The blood will be automatically absorbed.

Confirmation Window

This is where you confirm if enough blood has been applied to the absorbent hole in the strip.

Test Strip Handle

Hold this part to insert the test strip into the slot.

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

ATTENTION:

The front side of test strip should face up when inserting test strip.

Test results might be wrong if the contact bar is not fully inserted into the test slot.


NOTE:

The **FORA D40g Plus** meter should only be used with **FORA D40** Test Strips. Using other test strips with this meter may produce inaccurate results.

SETTING THE METER


This device is equipped with voice function to guide you through the process of blood glucose and blood pressure testing.

Before using your meter for the first time or if you change the meter battery, check and update these settings. Make sure you complete the steps below and have your desired settings saved.

Entering the Setting Mode

Start with the meter off (no test strip inserted). Press **S** until the meter turns on.


1. Setting the date and the time

With the year flashing, press **M** until the correct year appears. Press **S**.

With the month flashing, press **M** until the correct month appears. Press **S**.


With the day flashing, press **M** until the correct day appears. Press **S**.

With the hour flashing, press **M** until the correct hour appears. Press **S**.


With the minute flashing, press **M** until the correct minute appears. Press **S**.


2. Selecting the speaking volume

There are seven (7) speaking volume options to choose from. Press **M** until the desired speaking volume appears. To confirm your selection, press **S**.

3. Choose a language

Press **M** to select L1 / L2. The default language for the meter is L1, which is English; L2 is Spanish. To confirm your selection, press **S**.

4. Setting the auto-send

Press **M** to select the auto-send On or OFF. Press **S**. This function is for the cellular connectivity. If "On" is selected, your result will be transmitted right after each test.

5. Deleting the memory

If you do not wish to delete the saved results, press **S** to skip this step.

Otherwise, as the "dEL" and the "**M**" symbol are flashing on the display, and all of the user numbers are also flashing, press **M** twice to delete ALL results.

Press **S** to turn off the meter.

The voice guide repeats all settings before shutdown.

NOTE:


- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

THE FOUR MEASURING MODES


The meter provides you with four modes for measuring, General, AC, PC and QC.

MODES	USE WHEN
General Tests (not displayed)	any time of day without regard to time since last meal
AC	no food intake for at least 8 hours
PC	2 hours after a meal
QC	testing with the control solution

You can switch between each mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter.

2. Press **M** to switch between General, AC, PC and QC mode.

Control Solution Testing

FORA Control Solution contains a known amount of glucose that reacts with test strips and is used to ensure your meter and test strips are working together correctly.


Do a control solution test when:

- you first receive the meter,
- at least once a week to routinely check the meter and test strips,
- you begin using a new vial of test strips,
- you suspect the meter or test strips are not working properly,
- your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
- you want to practice the testing process; or
- you have dropped or think you may have damaged the meter.

Performing a Control Solution Test

1. Insert the test strip to turn on the meter


Insert the test strip into the meter. Wait for the meter to display “” and “”.

2. Press **M** to mark this test as a control solution test


WARNING:

When doing the control solution test, you have to mark it so that the test result will **NOT** mix with the blood glucose **TEST RESULTS** stored in the memory. Failure to do so will mix up the blood glucose test results with the control solution test results in memory.

3. Apply control solution

Shake the FORA control solution vial thoroughly before use. Squeeze out first drop and wipe it off, then squeeze out another drop and place it on the tip of the vial cap. Hold the meter to move the absorbent hole of the test strip to touch the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto a strip.

4. Read and compare the result

After counting down to 0, the control solution test result will appear on the display. Compare this result with the range printed on the test strip vial and it should fall within this range. If not, please read the instructions again and repeat the control solution test.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact the local customer service or place of purchase for help.

NOTE:

- The control solution range printed on the test strip vial is for control solution use only. It is not a recommended range for your blood glucose level.
- See the **MAINTENANCE** section for important information about your control solutions.
- It is recommended that the control solution test should be done at room temperature (68 °F to 77 °F/20 °C to 25 °C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.

TESTING WITH BLOOD SAMPLE

WARNING:

To reduce the chance of infection:

- Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets or the lancing device.
- Wash and dry your hands thoroughly before and after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

Sharing the lancing device and lancets may increase the risk of spreading infectious diseases. Lancing device must not be used for more than one person.

Preparing the Lancing Device for Blood Testing

Please follow the instructions in the lancing device insert for collecting a blood sample.

Preparing the Puncture Site

Stimulating blood perfusion by rubbing the puncture site before blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site was rubbed prior to blood extraction, the difference was significantly reduced.


Please follow the suggestions below before obtaining a drop of blood:

- Wash and dry your hands before starting.
- Select the puncture site either at fingertips or another body parts (please see section “Alternative Site Testing” (AST) on how to select the appropriate sites).
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**
- Rub the puncture site for about 20 seconds before penetration.
- Use a clear cap while setting up the lancing device.

• **Fingertip testing**

Press the lancing device's tip firmly against the lower side of your fingertip. Press the release button to prick your finger, then a click indicates that the puncture is complete.

- **Blood from sites other than the fingertip**

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm or hand, avoid lancing the areas with obvious veins because of excessive bleeding.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Please consult your health care professional before you begin AST.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Performing a Blood Glucose Test

1. Insert the test strip to turn on the meter



Wait for the meter to display “” and “”.

2. Select the appropriate measuring mode by pressing **M**

For selecting the measurement mode, please refer to the “FOUR MEASURING MODES”.

3. Obtaining a blood sample

Use the pre-set lancing device to puncture your desired site. After penetration, discard the first drop of blood with a




clean tissue or cotton. Gently squeeze the punctured area to obtain another drop of blood. Be careful **NOT** to smear the blood sample.

The volume of blood sample must be at least **0.5** microliter of volume. (● actual size)


4. Apply the sample

Hold the blood drop to touch the absorbent hole of the test strip. Blood will be drawn in and after the confirmation window is completely filled, the meter begins counting down.

NOTE:

- Do not press the punctured site against the test strip or try to smear the blood.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to start a new test.
- The confirmation window should be filled with blood before the meter begins to count down. **NEVER** try to add more blood to the test strip after the drop of blood has moved away. **Discard the used test strip and retest with a new one.**
- If you have trouble filling the confirmation window, please contact your health care professional or the local customer service for assistance.

5. Read your result

The result of your blood glucose test will appear after the meter counts down to 0. This blood glucose result will automatically be stored in the memory.

6. Eject the used test strip and remove the lancet

To eject the test strip, point the strip at a disposal container for sharp objects. The meter will switch itself off automatically after the test strip is ejected.

Always follow the instructions in the lancing device insert when removing the lancet.


WARNING:

- The used lancet and test strip may be potentially biohazardous. Please discard them carefully.
- Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination. For more information, please refer to the “Cleaning and Disinfection” section.

Alternative Site Testing

IMPORTANT: There are limitations with AST (Alternative Site Testing). Please consult your health care professional before you perform AST.

- Alternative site testing should not be used for calibrating continuous glucose meters (CGMs).
- Alternative site testing should not be used for insulin dose calculations.

What is AST?

Alternative site testing (AST) means that people use parts of the body other than the fingertips to check their blood glucose levels. This system allows you to test on the palm, the forearm, the upper arm, the calf or the thigh with the equivalent results to fingertip testing.

What is the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, since nerve endings are not so condensed, you will not feel as much pain at the fingertips.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after meal, physical exercise, or any other event, **take the blood sample from your finger only.**

We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or more than 2 hours after the last meal.
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not match the way you feel.
- Your routine glucose results often fluctuate.

TESTING YOUR BLOOD PRESSURE

The meter provides you two different ways to measure your blood pressure. Select from the options below:

- Single measurement

Perform an individual blood pressure measurement.

- Average measurement


Automatically perform three (3) consecutive blood pressure measurements and display the final average result.

Before Measurement


- Avoid caffeine, tea, alcohol and tobacco at least 30 minutes before measurement.
- Wait 30 minutes after exercising or bathing before measurement.
- Sit or lie down for at least 10 minutes before measurement.
- Do not measure when feeling anxious or tense.
- Take a 5-10 minute break between measurements. This break can be longer if necessary, depending on your physical conditions.
- Keep the records for your doctor as reference.
- Blood pressure varies between each arm. Always measure your blood pressure on the same arm.

- To take a blood pressure measurement after performing a blood glucose test, make sure that the test strip has been removed from the meter.


Fitting the Cuff Correctly

1. Connect the air plug of the tubing to the air jack of the meter.


2. Assemble the cuff as shown on the right. The smooth surface should be inside the cuff loop and the metal D-ring should not touch your skin.

3. Stretch your left (or right) arm in front of you with your palm facing up. Slide and place the cuff onto your arm with the air tube and artery mark region (in red) toward the lower arm.

Wrap and tighten the cuff above your elbow. The red line on the edge of the cuff should be approximately **0.8 to 1.2 inches (2 to 3 cm)** above your elbow. Align the tube over the main arteries on the inside.

4. Leave a little free space between the arm and the cuff, you should be able to fit 2 fingers between them. Clothing must not restrict the arm. Remove all clothing covering or constricting the measurement arm.

- The range index of cuff should fall into this range.

5. Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

Proper Measurement Position

1. Sit down for at least 10 minutes before measuring.

2. Place your elbow on a flat surface. Relax your hand with the palm facing up.

3. Make sure the cuff is about the same height as the location of your heart.

Remain still and do not talk or move during the measurement.

4. Measurement is in progress.


After the meter is turned on, the cuff will begin to inflate automatically.

Taking Measurements

This meter provides you two different ways to measure your blood pressure.

• Taking a Single Measurement

Always apply the pressure cuff before turning on the meter.

1. Press . All the LCD symbols will appear. Then the cuff will begin to inflate automatically.

2. The heart symbol “♥” will flash when a pulse is detected during the inflation.

3. After the measurement, the meter displays the systolic pressure, diastolic pressure and pulse rate.


4. Press  to switch off, or it will switch off automatically if left idle for 3 minutes.

NOTE:

- If you press  during measurement, the meter will be turned off.
- If the pulse rate symbol is shown as “” instead of “♥”, the appearance of symbol means the heart beat period deviating from the mean heart beat by greater than 25%.

• Measurement Average Mode

Always apply the pressure cuff before turning on the meter.

1. Press . The meter will turn on and enter the averaging mode. Then the cuff will begin to inflate automatically.

2. After the first measurement is complete, the meter will start counting down before the second measurement begins. The number on the right represents the remaining countdown between each measurement. The meter will take three (3) measurements consecutively with an interval of 20 seconds.


3. After taking three measurements, the results are averaged to produce the blood pressure measurement. Press  to turn off the meter.


METER MEMORY

Your meter stores the 864 most recent test results along with respective dates and times in its memory. To recall the memory, **start with the meter off.**



Reviewing Test Results

1. Press and release **M.** " **M** " will appear on the display.


2. Press **M.** Review all test results stored in the meter. Results display according to the time the measurement was taken. The latest measurement displays first. Press **M** to view older results, and press **S** to view newer results.

3. Exit the memory. Press **⏻** and the meter will turn off.

Reviewing Day Average Results

1. Press and hold **M for 3 seconds** until the “**DAY AVERAGE**” appears. Release **M** and then your blood glucose 7-day average result will appear on the display.

2. Press **M.** Review your 14-, 21-, 28-, 60- and 90- day average.

3. Press **S.** Switch to review your blood pressure 7-, 14-, 21-, 28-, 60- and 90-day average.

4. Exit the memory. Press **⏻** and the meter will turn off.

NOTE:

- Any time you wish to exit the memory, press **⏻** or leave it without any action for 3 minutes. The meter will turn off automatically.
- Control solution results are **NOT** included in the day average.
- If using the meter for the first time, “- - -” displays when you recall the test results or review the average result. This indicates that there is no test result in the memory.
- The blood glucose test results of “L” and “H” will **NOT** be calculated in the day average results.


DOWNLOADING RESULTS

Data transmission via cellular connectivity

The meter has completed the factory setting and pre-installed the SIM card before ex-factory. Do **NOT** remove the SIM card from the meter by yourself.

1. How to initiate cellular connectivity

- (a) with the meter turned off, press  to initiate cellular connectivity of batch upload for current user with test results which have not been uploaded.
- (b) setting the auto-send On, the meter will initiate cellular connectivity after measurement to send the current new data,
- (c) in memory mode, press  to initiate cellular connectivity to upload data shows on the screen.

2. During transmission

With the meter turned off, press  to initiate the data transmission from the meter. "PCL" will appear on the display first and the RF indicator will light up in blue. Then the test result under transmission will show on the screen.

When the RF indicator is blinking in blue, the data transmission has started. When the RF indicator shows an orange light, the cellular connectivity service is online.



Press  and the meter will turn off.

NOTE:

- If the SIM card is not inserted properly, the “PCL” will start blinking, and the meter will turn off automatically after blinking for 30 seconds.
- If the meter displays a message like the one on the right while connecting to the server, please try again. If the problem still persists, contact the place of purchase for help.
- While the meter is doing a data transmission, it is unable to perform tests.

Signal Quality Check



With the meter switched off, press and hold  for 10 seconds to begin the signal quality check. Wait until 'PCL' appears on the screen.

The number displayed in the middle of the screen represents the signal quality:

44 to 94: Acceptable

95 to 140: Weak Signal

Poor signal quality may be caused by:

- the service is busy; or
- there is no base station nearby.

Poor signal quality may affect the data transfer process. If the signal quality is poor, move the meter to a spot with better reception and try again. Make sure the meter is not close to or nearby a microwave, mobile devices, or any other devices that may interrupt the signal.

When signal quality check is complete, press  to turn off the meter.


MAINTENANCE

Battery

Your meter comes with one 3.7V lithium-ion rechargeable battery.

Low Battery Signal

The meter will display the two messages below to alert you when the meter power is getting low.

1. The “” appears along with display messages: The meter is functional and the result remains accurate, but it is time to charge the battery.

If the power is not enough to do a test, the symbol starts blinking. Please charge the Li-ion rechargeable battery immediately.

2. The power indicator lights up in red.


Replacing the Battery

All rechargeable batteries are consumables and have a limited lifespan. As the Li-ion battery chemically ages, the amount of charge it can hold diminishes, resulting in shorter amounts of time before a device needs to be recharged. Replace a new Li-ion battery immediately when the battery's amount of charge diminishes or the battery expands. Please contact the supplier to purchase Li-ion battery separately.

To replace the Li-ion battery, make sure that the meter is turned off.

WARNING:

- FORA D40g Plus can use with one Li-ion rechargeable battery **ONLY** (battery pack no.S14500 1S4P).
- Use only new Li-ion rechargeable battery of the specified size and type. Wrong types of battery may cause explosion.
- The AC Power Adaptor is only for charging the Li-ion rechargeable battery. The AC Power Adaptor alone can't provide sufficient power for machine operation.

- 1.** Press the edge of the battery cover and lift it up to remove.
- 2.** Remove the old Li-ion battery and replace with new one 3.7V rechargeable Li-ion battery.
- 3.** Close the battery cover. If the Li-ion battery is inserted correctly, you will hear a “beep” afterwards.

CAUTION: Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.

The Li-ion battery is about the same shape and size as the four AA alkaline batteries. If you are installing a Li-ion battery, make sure the battery is correctly inserted in the compartment.

NOTE:

- Replacing the batteries does not affect the test results stored in memory.
- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Batteries may leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e. 3 months or more).
- Properly dispose of the batteries.

Recharging the Battery

The power adapter can be used as the battery charger when you need to recharge Li-ion battery. To recharge the 3.7V Li-ion battery, connect AC adapter to the meter and a power source.

The power indicator will light up in green, which indicates the battery is undergoing a recharge. It should take approximately 5 hours to be fully charged.

If the power indicator is blinking in red when connecting AC adapter to a power source, this indicates the meter is in a pre-charging stage. If the power indicator starts blinking in green, this indicates the recharge is nearly finished. The power indicator will diminish when the recharge completes.

CAUTION:

- The Li-ion battery must be recharged with the **FORA D40g Plus** meter by connecting AC adapter to an electrical outlet. Only this 3.7V rechargeable Li-ion battery (battery pack no. S14500 1S4P) can be recharged by **FORA D40g Plus** with the connection of AC adapter. Any other kind of rechargeable battery is not allowed or it may damage the meter.

Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash hands thoroughly with soap and water before and after use.

Why the cleaning and disinfection should be performed

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter

All surface of meter if visibly soiled must be physically cleaned to remove gross soil. Disinfect the meter between each patient to prevent infection.

How to clean and disinfect the meter



The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt, blood, or any other body fluid with the wipe. Use a second wipe to disinfect the meter by following the disinfecting procedure below. Do **NOT** use organic solvents to clean the meter.

We recommend for meter cleaning and disinfection you should use the disinfecting wipes/towelettes from below. It has been shown to be safe for use with the **FORA D40g Plus** Blood Glucose plus Blood Pressure Monitoring System.

▶ Micro-Kill+™ (Micro-Kill Plus™) by Medline (EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com

Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.
2. Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep meter wet with disinfection solution for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe.
Do not use the wipe to clean or disinfect the pressure cuff.
3. Remove the wipe. Allow the meter surface to dry completely.
4. Discard the used wipes and never reuse them.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact local customer service at 1-888-307-8188 for assistance.

Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination. Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact local customer service at 1-888-307-8188 for assistance.

This device has been validated to withstand up to 10,000 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 5 cleaning and

disinfection cycles per day over 5 years, the expected life of the device. The meter should be replaced after the validated number of cleaning and disinfection cycles or the warranty period, whichever comes first.

Stop using the meter if you see any signs of deterioration. For example:

- meter can not be turned on,
- LCD display cracks or becomes cloudy,
- buttons no longer function,
- meter outer casing cracks,
- data can not be transmitted to PC,
- color or paint/printing on housing is abnormal, or
- scratches or abrasions on meter are higher than acceptable.

Please contact customer service for a replacement meter if any of the signs of deterioration are noticed.

NOTE:

- Do **NOT** clean and disinfect the meter while performing tests.
- Do **NOT** flush your meter with water.
- Do **NOT** wash or iron the pressure cuff.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot, battery compartment, or strip-ejection button.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

Meter Disposal

The used meter should be treated as a contaminant that may carry a risk of infection during measurement. The batteries in this used meter should be removed and the meter should be disposed properly.

Meter Storage

- Storage condition: -4°F - 140°F (-20°C - 60°C), between 10% to 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

Caring for Your Test Strips

- Storage condition: 39.2°F - 104°F (4°C - 40°C), between 10% to 85% relative humidity. Do not freeze.
- Store your test strips in their original vial only. Do not transfer to other container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on the strip vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiration date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

Important Control Solution Information

- Use only FORA control solutions with your meter.
- Do not use the control solution beyond the expiration date or 3 months after first opening. Write the opening date on the control solution vial and discard the remaining solution after 3 months.
- It is recommended that the control solution test should be done at room temperature (68°F - 77°F/20°C - 25°C). Make sure your control solution, meter, and test strips are at this specified temperature range before testing.
- Shake the vial before use, discard the first drop of control solution, and wipe off the dispenser tip to ensure a pure sample and an accurate result.
- Store the control solution tightly closed at temperatures between 36°F and 86°F (2°C to 30°C). Do NOT freeze.

Caring for Your Pressure Cuff

- To clean the cuff, wipe it with a soft, moistened cloth and soap.
- Do not wash the cuff.
- Do not immerse the cuff in water or other liquids.
- Do not use petrol, thinners or similar solvents to clean the cuff.

Calibration and Service

The accuracy of this blood pressure meter has been carefully tested and is designed for a long service life. It does not require calibration. Please consult your local authorized FORA distributor or dealer for any question.

SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. Do not attempt to repair by yourself and never try to disassemble the meter under any circumstances.

Result Readings

MESSAGE	WHAT IT MEANS
BLOOD GLUCOSE	
 The image shows a digital display on a blood glucose meter. At the top, it displays '10-12' on the left and '10:00' on the right. In the center, the word 'Lo' is displayed in a large, bold font. There are small icons above the 'Lo' and '10:00'.	Appears when your result is below measurement limit, which is less than 20 mg/dL (1.1 mmol/L).
 The image shows a digital display on a blood glucose meter. At the top, it displays '10-12' on the left and '10:00' on the right. In the center, the number '250' is displayed in a large, bold font. There are small icons above the '250' and '10:00'.	Appears when your result is equal to or higher than 240 mg/dL (13.3 mmol/L). This indicates the possibility of ketone accumulation for type 1 diabetes. Please seek medical assistance immediately.
 The image shows a digital display on a blood glucose meter. At the top, it displays '10-12' on the left and '10:00' on the right. In the center, the letters 'Hi' are displayed in a large, bold font. There are small icons above the 'Hi' and '10:00'.	Appears when your result is higher than the limit of measurement, which is higher than 600 mg/dL (33.3mmol/L).

Error Messages

MESSAGE	WHAT IT MEANS	WHAT TO DO
Err 1	Problem with the cuff or meter.	Please contact local customer service for help.
Err 2~Err 8	Problem with the measurement.	Refit cuff tightly and correctly. Relax for 5 minutes and repeat the measurement. If the error still remains, please contact the local customer service for help.
Err 10~Err 20	Problem with the meter.	Please contact the local customer service for assistance.
E-U	Strip has been used.	Repeat the test with a new strip.
E-F	Test strip is removed while counting down, or insufficient blood volume.	Review the instructions and repeat the test with a new strip. If problem persists, contact the local customer service for assistance.
E-t	Environmental temperature is outside the system's operational range.	System operational range is 50°F to 104°F (10°C to 40°C). Repeat the test after the meter and test strip have reached the above temperature.
	Battery is low.	Recharge the battery immediately.

Error Messages of Cellular Connectivity

MESSAGE	WHAT IT MEANS	WHAT TO DO
PCL	No SIM card inserted.	Please contact customer service for assistance.
PCL Err 27–30 PCL Err 40–45 PCL Err 50 PCL Err 62–69	Problem with upload.	Step 1 : Please start over and try again. Step 2 : If problem still persists, please contact local customer service for assistance.
PCL Err 38 PCL Err 72-76	Problem with setting.	Please contact customer service for assistance.
PCL Err 39	Problem with hardware.	Turn off the device and start over. Follow the instructions and try again.
PCL Err 2-26 PCL Err 32-37 PCL Err 46–61 PCL Err 84–99 PCL Err 119–127	Problem with network.	Step 1 : Please start over and try again. Step 2 : If problem still persists, please contact local customer service for assistance.
PCL Err 70–71	Problem with hardware.	Step 1 : Remove the battery. Step 2 : Install the battery. Step 3 : Recharge the battery.
PCL Err77	Problem with network.	Step 1 : Turn off the device and reconnect to the service. Step 2 : Change to a spot with a better quality of connection and try again. Step 3 : If problem still persists, please contact local customer service for assistance.

Troubleshooting

Blood Glucose Measurement

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer services.

2. If the test does not start after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer services.
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood sample.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature 68° F to 77° F (20° C to 25° C) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer services.

Blood Pressure Measurement

1. If no display appears after pressing .

POSSIBLE CAUSE	WHAT TO DO
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

2. If the heart rate is higher/lower than user's average:

POSSIBLE CAUSE	WHAT TO DO
Movement during measurement.	Repeat measurement and remain still during test.
Measurement taken just after exercise.	Rest at least 30 minutes before repeating measurement.

3. If the result is higher/lower than user's average measurement:

POSSIBLE CAUSE	WHAT TO DO
May not be in correct position while measuring.	Adjust to the correct position to measure.
Blood pressure naturally varies from time to time.	Keep in mind for next measurement.

DETAILED INFORMATION

Reference Values

● Blood Glucose

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association. Classification and Diagnosis of Diabetes: Standards of Medical Care in Diabetes—2022 Jan; 45(Supplement 1): S17-S38.

Please consult your doctor to determine a target range that works best for you.

● Blood Pressure

Clinical studies show that adult diabetes is often accompanied by elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment*¹.

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction.

The recommended blood pressure range is as below:

Classification	Systolic Pressure (mmHg)	Diastolic Pressure (mmHg)
Hypotension* ²	< 90	< 60
Normal* ³	< 120	< 80
Elevated* ³	120 -129	< 80
Stage 1 Hypertension* ³	130 -139	80 -89
Stage 2 Hypertension* ³	≥ 140	≥ 90
Hypertensive Crisis* ³	>180	>120

***1:** National Heart, Lung, and Blood Institute, Diseases and Conditions.

***2:** The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure. NIH Publication. 2003. No. 03-5233.

***3:** New ACC (American College of Cardiology) and AHA (American Heart Association) guidelines for the detection, prevention, management and treatment of high blood pressure (2017).

SPECIFICATIONS

- **System Performance**

Power Source:

One 3.7V lithium-ion rechargeable battery (battery pack no. S14500 1S4P)

Size of Meter w/o Cuff:

147 mm (L) x 105 mm (W) x 80 mm (H),

Weight: 460g approx. (with batteries)

Memory: Maximum 864 memory records

Power Saving:

Automatic power off if system is idle for 3 minutes (normal mode) or 5 minutes (RF mode)

System Operating Condition:

50°F to 104°F (10°C to 40°C), 10 - 93% RH

Meter Storage Condition:

-4°F to 140°F (-20°C to 60°C), 10 - 95% RH

Strip Storage Condition:

39.2°F to 104°F (4°C to 40°C), 10 - 85% RH

Power Supply Input: DC +6V / 1A (max) via Power Plug

Expected Service Life: 3 years

● **Blood Glucose Measurement Performance**

Measurement Unit: mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

Precision: $\pm 5\%$ (CV)

Accuracy:

± 15 mg/dL (0.83 mmol/L) when glucose ≤ 75 mg/dL (4.16 mmol/L)

$\pm 20\%$ when glucose > 75 mg/dL (4.16 mmol/L)

Ketone Warning: ≥ 240 mg/dL (13.3 mmol/L)

● **Blood Pressure Measurement Performance**

Measurement Unit: mmHg

Pressure Range: 0 - 300 mmHg

Systolic Measurement Range: 60 to 250 mmHg

Diastolic Measurement Range: 30 to 180 mmHg

Pulse Rate Measurement Range: 30 to 199 beats / minute

Maximum Inflation Pressure: 280 mmHg

Accuracy of Pressure: ± 3 mmHg

Accuracy of Pulse Rate: $\pm 4\%$ of reading

This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC/EN 61010-1, IEC/EN 61010-2-101, IEC/EN 61326-2-6.

Reference to Standards:

- IEC 60601-1 General requirement for safety
- IEC 60601-1-2 Requirements for EMC
- AAMI/ANSI /IEC 80601-2-30, ANSI/AAMI/ISO 81060-2, NIBP requirements

FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.

15.105(b)

Federal Communications Commission (FCC) Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including interference that may cause undesired operation of the device.

FCC RF Radiation Exposure Statement:

1. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
2. This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

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Distributed by Foracare, Inc. / Distribuido por Foracare, Inc.