

BINSON'S

Medical Equipment & Supplies

Since 1953



G20

BLOOD GLUCOSE MONITORING SYSTEM

FEATURES & BENEFITS

- Strip Ejector
- No coding required
- Tiny Blood Sample (0.7 µL)
- 448 Test Memory

SPECIFICATIONS

Blood Test Volume: 0.7 µL

Reaction Time: 7 seconds

Measurement Range: 20 ~ 600 mg/dL
(1.1~33.3 mmol/L)

Memory Capacity: 448 results with date & time

Ketone Warning: Yes \geq 240mg/dL

Day Average: 7/14/21/28/60/90 Days

Power Source: One Lithium Battery (CR2032)

Data Output: Standard RS232 PC Interface



Attention Medicare Beneficiaries! You can receive your FORA G20 meter and supplies now from Binson's!

HOME DELIVERY DIABETES SUPPLIER

Call 1-855-307-2973 to place your order. Make the change to Binson's today!