

# BINSON'S

Medical Equipment & Supplies

Since 1953



## G20

BLOOD GLUCOSE MONITORING SYSTEM

## FEATURES & BENEFITS

- Strip Ejector
- No coding required
- Tiny Blood Sample (0.7 µL)
- 448 Test Memory

## SPECIFICATIONS

**Blood Test Volume:** 0.7 µL

**Reaction Time:** 7 seconds

**Measurement Range:** 20 ~ 600 mg/dL  
(1.1~33.3 mmol/L)

**Memory Capacity:** 448 results with date & time

**Ketone Warning:** Yes  $\geq$  240mg/dL

**Day Average:** 7/14/21/28/60/90 Days

**Power Source:** One Lithium Battery (CR2032)

**Data Output:** Standard RS232 PC Interface



**Attention Medicare Beneficiaries! You can receive your FORA G20 meter and supplies now from Binson's!**

## HOME DELIVERY DIABETES SUPPLIER

Call 1-855-307-2973 to place your order. Make the change to Binson's today!