

BINSON'S

Medical Equipment & Supplies

Since 1953

FORA® V10

BLOOD GLUCOSE MONITORING SYSTEM



FEATURES & BENEFITS

- Friendly talking function (English/Spanish)
- No coding required
- Alternative site testing approved
- Data transmission available for better management

SPECIFICATIONS

Blood Test Volume: 0.7 µL
Reaction Time: 7 seconds
Measurement Range: 20 ~ 600 mg/dL (1.1~33.3 mmol/L)
Memory Capacity: 450 results with date & time
Ketone Warning: Yes ≥ 240 mg/dL
Day Average: 7/14/21/28/60/90 Days
Power Source: Two AAA Alkaline Battery
Data Output: Standard RS232 PC Interface



Attention Medicare Beneficiaries! You can receive your FORA V10 meter and supplies now from Binson's!

HOME DELIVERY DIABETES SUPPLIER

Call 1-855-307-2973 to place your order. Make the change to Binson's today!